

Section B

Positioning Accessories

- B 1.- Positioning Support Pads and Hardware (pg. 30)
 - B 1.1.- Positioning Pads (pg. 30)
 - B 1.1.1.- Lateral Trunk Support Pads (pg. 31)
 - B 1.1.1.1.- Flat Lateral Trunk Support Pads (pg. 31)
 - B 1.1.1.2.- Contoured Lateral Trunk Support Pads (pg. 32)
 - B 1.1.2.- Lateral Hip/Thigh Support Pads (pg. 33)
 - B 1.1.2.1.- Flat Lateral Hip/Thigh Support Pads Mounted From The Back (pg. 33)
 - B 1.1.2.2.- Contoured Lateral Hip/Thigh Support Pads Mounted From The Back (pg. 34)
 - B 1.1.2.3.- Flat Lateral Hip/Thigh Support Pads Mounted From The Seat (pg. 35)
 - B 1.1.2.4.- Contoured Lateral Hip/Thigh Support Pads Mounted From The Seat (pg. 36)
 - B 1.1.2.5.- Flat Lateral Hip/Thigh Support Pads Mounted From Armrests (pg. 37)
 - B 1.1.3.- Lateral Knee Support Pads (pg. 38)
 - B 1.1.3.1.- Flat Lateral Knee Support Pads (pg. 38)
 - B 1.1.3.2.- L-Shaped Lateral Knee Support Pads (pg. 39)
 - B 1.1.3.3.- Round Shaped Lateral Knee Support Pads (pg. 40)
 - B 1.1.4.- Medial Knee Support Pads (pg. 41)
 - B 1.1.4.1.- Medial Knee Support Pads, Wedge Shaped (pg. 41)
 - B 1.1.4.2.- Medial Knee Support Pads, L-Shaped Wedge (pg. 42)
 - B 1.1.4.3.- Medial Knee Support Pads, Oval Shaped (pg. 42)
 - B 1.1.4.4.- Medial Knee Support Pads, L-Shaped Oval (pg. 43)
 - B 1.2.- Mounting Hardware for Positioning Pads (pg. 44)
 - B 1.2.1.- Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports, Lateral Knee Supports and Medial Knee Supports (pg. 44)
 - B 1.2.1.1.- 90 Degree Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports (pg. 44)
 - B 1.2.1.2.- Offset Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports (pg. 46)
 - B 1.2.1.3.- Removable Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports (pg. 46)
 - B 1.2.1.4.- Single Knob Assembly for Mounting Hardware of Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports (pg. 47)
 - B 1.2.1.5.- Rubber Coated Metal Straps for Lateral Hip/Thigh Supports, Lateral Knee Supports (pg. 48)
 - B 1.2.1.6.- Swing Away Hardware for Lateral Trunk Supports (pg. 48)
 - B 1.2.1.7.- Hardware for Medial Knee Supports (pg. 49)
 - B 1.2.1.8.- Hardware for Round Lateral Knee Supports (pg. 49)
 - B 1.2.2.- Pre-Configured Mounting Hardware for Positioning Pads (pg. 50)
 - B 1.2.2.1.- Pre-Configured Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports, Lateral Knee Supports and Medial Knee Supports (pg. 50)
 - B 1.2.2.2.- Pre-Configured Mounting Hardware for Lateral Hip/Thigh Supports secured from The Armrest (pg. 53)
 - B 1.2.2.3.- Pre-Configured Mounting Hardware for Round Lateral Knee Supports (pg. 54)
 - B 1.2.2.4.- Pre-Configured Mounting Hardware for Medial Knee Supports (pg. 54)
- B 2.- Positioning Belts, Straps and Pads (pg. 55)
 - B 2.1.- Anterior Trunk Supports (pg. 55)
 - B 2.1.1.- BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support (pg. 56)
 - B 2.1.2.- BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support (pg. 57)
 - B 2.1.3.- BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure (pg. 57)
 - B 2.1.4.- BioForm Stretch (Black Neoprene) "H"-Style Anterior Trunk Support (pg. 58)
 - B 2.1.5.- BioForm Stretch (Black Neoprene) 6 Point Trunk Support (pg. 58)
 - B 2.1.6.- BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support (pg. 59)
 - B 2.1.7.- BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support (pg. 59)
 - B 2.1.8.- BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure (pg. 60)

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

- B 2.1.9.- BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support (pg. 60)
- B 2.1.10.- BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support (pg. 61)
- B 2.1.11.- BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support (pg. 61)
- B 2.1.12.- BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support (pg. 62)
- B 2.1.13.- BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure (pg. 62)
- B 2.1.14.- BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support (pg. 63)
- B 2.1.15.- BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support (pg. 63)
- B 2.1.16.- BioForm Black Ballistic Nylon Standard Anterior Trunk Support (pg. 64)
- B 2.2.- BioForm Abdominal Binders (pg. 64)
- B 2.3.- BioForm Positioning Belts (pg. 65)
 - B 2.3.1.- BioForm Positioning Belts, Velcro Closure (pg. 65)
 - B 2.3.2.- BioForm Positioning Belts, Fastex Buckle Closure (pg. 65)
 - B 2.3.3.- BioForm Positioning Belts, Auto Style Buckle Closure (pg. 65)
- B 2.4.- BioForm Positioning Belts with Split Pads (pg. 66)
 - B 2.4.1.- BioForm Positioning Belts, Split Pads, 2 each, Fastex Buckle Closure (pg. 66)
 - B 2.4.2.- BioForm Positioning Belts, Split Pads, 2 each, Auto-Style Buckle Closure (pg. 66)
 - B 2.4.3.- BioForm 4-Points Positioning Belts, Split Pads, 2 each, Auto-Style Buckle Closure (pg. 66)
- B 2.5.- BioForm Chest Positioning Belts (pg. 67)
 - B 2.5.1.- BioForm Non-Stretch Chest Positioning Belts (pg. 67)
 - B 2.5.2.- BioForm Stretch (Black Neoprene) Chest Positioning Belts (pg. 67)
- B 2.6.- BioForm Calf Straps (pg. 68)
- B 2.7.- BioForm Thigh Straps (pg. 68)
- B 2.8.- BioForm Padded Sleeves (pg. 69)
 - B 2.8.1.- BioForm Black Padded Sleeves For 2” Belts (pg. 69)
 - B 2.8.2.- BioForm Black Padded Sleeves For 1½” Belts (pg. 69)
- B 3.- Lower Extremity Supports (pg. 70)**
 - B 3.1.- BioForm Foot Positioners (pg. 70)
 - B 3.1.1.- BioForm Ankle Positioners (pg. 70)
 - B 3.1.2.- BioForm Multi-Slots Footplate Extenders (pg. 71)
 - B 3.1.3.- BioForm Ankle Positioner with Adjustable Plate (pg. 71)
 - B 3.1.4.- BioForm Foot Positioner (pg. 71)
 - B 3.2.- BioForm Footplates (pg. 72)
 - B 3.2.1.- BioForm One Piece, Footplates (pg. 72)
 - B 3.2.1.1.- BioForm One Piece, Fixed Footplate (pg. 72)
 - B 3.2.1.2.- BioForm Padded One Piece, Fixed Footplate (pg. 72)
 - B 3.2.1.3.- BioForm One Piece, Flip-Up Footplate (pg. 73)
 - B 3.2.1.4.- BioForm Padded One Piece, Flip-Up Footplate (pg. 73)
 - B 3.2.2.- BioForm Footplates Extenders (pg. 74)
 - B 3.2.2.1.- BioForm Footplates Extender (pg. 74)
 - B 3.2.2.2.- BioForm Padded Footplates Extender (pg. 74)
 - B 3.2.3.- BioForm Footplate Pads (pg. 75)
 - B 3.2.3.1.- BioForm Soft Footplate Pads (pg. 75)
 - B 3.2.3.2.- BioForm Footplate Pads (pg. 75)
 - B 3.3.- BioForm Sandals and Straps (pg. 76)
 - B 3.3.1.- BioForm Sandals with one Strap w/Velcro (pg. 76)
 - B 3.3.2.- BioForm Sandals with one Strap w/Buckle (pg. 76)
 - B 3.3.3.- BioForm Sandals with one Leather Strap (pg. 77)
 - B 3.3.4.- BioForm Sandals with two Straps w/Velcro (pg. 77)
 - B 3.3.5.- BioForm Sandals with two Straps w/Buckle (pg. 77)
 - B 3.3.6.- BioForm Sandals with two Leather Straps (pg. 78)
 - B 3.4.- Solid Heel Cups (pg. 78)
 - B 3.4.1.- BioForm Solid Heel Cups with Straps w/ Velcro (pg. 78)
 - B 3.4.2.- BioForm Solid Heel Cups with Straps w/ Fastex Buckle (pg. 79)
 - B 3.5.- BioForm Foot Support Accessories (pg. 79)
 - B 3.5.1.- BioForm Scuba Heel Loops (pg. 79)
 - B 3.5.2.- BioForm Ankle or Toe Straps (pg. 80)
 - B 3.5.2.1.- BioForm Ankle or Toe Straps, Velcro Closure (pg. 80)
 - B 3.5.2.2.- BioForm Ankle or Toe Straps, Fastex Buckle Closure (pg. 80)
 - B 3.5.2.3.- BioForm Padded Ankle or Toe Straps, Fastex Buckle Closure (pg. 80)
 - B 3.5.2.4.- Pads for Ankle or Toe Straps (pg. 81)
 - B 3.5.2.5.- BioForm Leather Straps (pg. 81)

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

- B 3.6.- BioForm Footboxes (pg. 82)
 - B 3.6.1.- BioForm, Double Leg (1 piece), Non-padded Footboxes (pg. 82)
 - B 3.6.1.1.- BioForm, Double Leg (1 piece), Non-padded Footboxes, No Hardware (pg. 82)
 - B 3.6.1.2.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Nuts and Bolts to secure from footplates. (pg. 83)
 - B 3.6.1.3.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Hardware to secure from Hangers. (pg. 83)
 - B 3.6.1.4.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Hardware to secure from Seat Pan. (pg. 84)
 - B 3.6.2.- BioForm, Single Leg (2 pieces), Non-padded Footboxes (pg. 84)
 - B 3.6.2.1.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, No Hardware (pg. 84)
 - B 3.6.2.2.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Nuts and Bolts to secure from footplates. (pg. 85)
 - B 3.6.2.3.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Hardware to secure from Hangers. (pg. 85)
 - B 3.6.2.4.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Hardware to secure from Seat Pan. (pg. 86)
 - B 3.6.3.- BioForm, Double Leg (1 piece), Padded Footboxes (pg. 86)
 - B 3.6.3.1.- BioForm, Double Leg (1 piece), Padded Footboxes, No Hardware (pg. 86)
 - B 3.6.3.2.- BioForm, Double Leg (1 piece), Padded Footboxes, Nuts and Bolts to secure from footplates. (pg. 87)
 - B 3.6.3.3.- BioForm, Double Leg (1 piece), Padded Footboxes, Hardware to secure from Hangers. (pg. 87)
 - B 3.6.3.4.- BioForm, Double Leg (1 piece), Padded Footboxes, Hardware to secure from Seat Pan. (pg. 88)
 - B 3.6.4.- BioForm, Single Leg (2 pieces), Padded Footboxes (pg. 88)
 - B 3.6.4.1.- BioForm, Single Leg (2 pieces), Padded Footboxes, No Hardware (pg. 88)
 - B 3.6.4.2.- BioForm, Single Leg (2 pieces), Padded Footboxes, Nuts and Bolts to secure from footplates. (pg. 89)
 - B 3.6.4.3.- BioForm, Single Leg (2 pieces), Padded Footboxes, Hardware to secure from Hangers. (pg. 89)
 - B 3.6.4.4.- BioForm, Single Leg (2 pieces), Padded Footboxes, Hardware to secure from Seat Pan. (pg. 90)
- B 3.7.- BioForm Calf Support (pg. 91)
 - B 3.7.1.- BioForm Soft Calf Panels (pg. 91)
 - B 3.7.1.1.- BioForm Stretch Soft Black Calf Panels with Fastex Buckle Closure (pg. 91)
 - B 3.7.1.2.- BioForm Stretch Soft Black Calf Panels with Velcro Closure (4+2 straps) (pg. 92)
 - B 3.7.1.3.- BioForm Non-Stretch Soft Black Calf Panels with Fastex Buckle Closure (pg. 92)
 - B 3.7.1.4.- BioForm Non-Stretch Soft Black Calf Panels with Velcro Closure (4+2 straps) (pg. 93)
 - B 3.7.2.- BioForm Solid Calf Panels (pg. 93)
 - B 3.7.2.1.- BioForm Solid Calf Panel, Fixed Mounting Hardware (pg. 93)
 - B 3.7.2.2.- BioForm Solid Calf Panel, Removable Hardware (pg. 94)
- B 3.8.- BioForm Anterior Knee Blocks (pg. 95)
- B 4.- Head Support Systems (pg. 96)
 - B 4.1.- Head Support Pads (pg. 96)
 - B 4.1.1.- Flat Headrest Pad (pg. 97)
 - B 4.1.2.- BioCurve Headrest Pad (pg. 97)
 - B 4.1.3.- BioForm Cervical Control Headrest Pad (pg. 97)
 - B 4.1.4.- BioForm Tri-Pad Headrest with Angle Adjustable Wings (pg. 98)
 - B 4.1.5.- BioForm Tri-Pad Headrest with Fixed Angle Wings (pg. 98)
 - B 4.1.6.- BioCurve Headrest Package (pg. 99)
 - B 4.2.- Headrest Support Modifications (pg. 99)
 - B 4.3.- Headrest Support Mounting Hardware (pg. 100)
- B 5.- Upper Extremity Supports (pg. 101)
 - B 5.1.- Upper Extremity Support Surface (pg. 101)
 - B 5.1.1.- BioDezine Upper Extremity Support Surface, Package (pg. 102)
 - B 5.1.2.- BioForm Upper Extremity Support Surface, Package (pg. 103)
 - B 5.1.3.- Upper Extremity Support Surface Materials (pg. 104)
 - B 5.1.4.- Upper Extremity Supports Surface Modifications (pg. 105)
 - B 5.1.5.- Upper Extremity Support Surface Mounting Hardware (pg. 106)
- B 6.- Miscellaneous Accessories (pg. 107)
 - B 6.1.- BioForm Joystick Handles (pg. 107)
 - B 6.1.1.- BioForm Goal Post Joystick Handles (pg. 107)
 - B 6.1.2.- BioForm 'T' Joystick Handles (pg. 107)

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.- Positioning Support Pads and Hardware

The Positioning Pads offered in this section enable the proper seating of the client. The extensive selection of hardware available for these pads facilitates the proper placement of the pads. Some of the hardware offered allows for the easy removal of the pad while other choices allow for the pad to be swung away.

B 1.1.- Positioning Pads



“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



B 1.1.1.- Lateral Trunk Support Pads

B 1.1.1.1.- Flat Lateral Trunk Support Pads

The Flat Lateral Trunk Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The upper trunk will be supported and stabilized allowing for more efficient use of the upper and lower extremities.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with. We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Flat Lateral Trunk Support Pads

Model	Description	Size	Dimension	
			H	D
BLD-NB	Flat Lateral Trunk Support Pads, pair	X-Small	3	4
BLS-NB	Flat Lateral Trunk Support Pads, pair	Small	4	4
BLM-NB	Flat Lateral Trunk Support Pads, pair	Medium	4	5
BLL-NB	Flat Lateral Trunk Support Pads, pair	Large	4	6
BLO-NB	Flat Lateral Trunk Support Pads, pair	X-Large	5	6
BLX-NB	Custom Flat Lateral Trunk Support Pads	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- BLM-NB becomes BLM-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


B 1.1.1.2.- Contoured Lateral Trunk Support Pads

The Contoured Lateral Trunk Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The contouring of the Lateral Trunk Support Pad will be beneficial by distributing weight/pressure over a greater surface area when the pad is used for postural correction and significant pressures will exist between the pad and consumer. The upper trunk will be supported and stabilized allowing for more efficient use of the upper and lower extremities.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs



Contoured Lateral Trunk Support Pads

Model	Description	Size	Dimension	
			H	D
BND-NB	Contoured Lateral Trunk Support Pads, pair	X-Small	3	4
BNS-NB	Contoured Lateral Trunk Support Pads, pair	Small	4	4
BNM-NB	Contoured Lateral Trunk Support Pads, pair	Medium	4	5
BNL-NB	Contoured Lateral Trunk Support Pads, pair	Large	4	6
BNO-NB	Contoured Lateral Trunk Support Pads, pair	X-Large	5	7
BNX-NB	Custom Contoured Lateral Trunk Support Pads	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- BNS-NB becomes BNS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.1.2.- Lateral Hip/Thigh Support Pads


B 1.1.2.1.- Flat Lateral Hip/Thigh Support Pads Mounted From The Back

The Flat Lateral Hip/Thigh Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts as well as keeping the pelvis centered on the seat. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The pelvis will be supported and stabilized allowing for more efficient use of the upper and lower extremities. By choosing from the Hip/Thigh Support product numbers below, the pads are to be mounted to the Solid Back Insert.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Flat Lateral Hip/Thigh Support Pads Mounted From The Back

Model	Description	Size	Dimension	
			H	D
BHD-NB	Flat Lateral Hip/Thigh Support Pads, pair	X-Small	3	4
BHS-NB	Flat Lateral Hip/Thigh Support Pads, pair	Small	4	5
BHM-NB	Flat Lateral Hip/Thigh Support Pads, pair	Medium	4	6
BHL-NB	Flat Lateral Hip/Thigh Support Pads, pair	Large	4	8
BHO-NB	Flat Lateral Hip/Thigh Support Pads, pair	X-Large	5	8
BHX-NB	Custom Flat Lateral Hip/Thigh Support Pads	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- BHS-NB becomes BHS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


B 1.1.2.2.- Contoured Lateral Hip/Thigh Support Pads Mounted From The Back

The Contoured Lateral Hip/Thigh Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts as well as keeping the pelvis centered on the seat. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The contouring of the Lateral Hip/Thigh Support Pad will be beneficial by providing a subtle cue to the Trochanteric area to assist in maintaining a properly centered pelvis. The pelvis will be supported and stabilized allowing for more efficient use of the upper and lower extremities. By choosing from the Hip/Thigh Support product numbers below, the pads are to be mounted to the Solid Back Insert.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Contoured Lateral Hip/Thigh Support Pads Mounted From The Back

Model	Description	Size	Dimension	
			H	D
BID-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	X-Small	3	4
BIS-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	Small	4	5
BIM-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	Medium	4	6
BIL-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	Large	4	8
BIO-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	X-Large	5	8
BIX-NB	Custom Contoured Lateral Hip/Thigh Support Pads, contoured at rear	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- BIS-NB becomes BIS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.1.2.3.- Flat Lateral Hip/Thigh Support Pads Mounted From The Seat

The Flat Lateral Hip/Thigh Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts as well as keeping the pelvis centered on the seat. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The pelvis will be supported and stabilized allowing for more efficient use of the upper and lower extremities. By choosing from the Hip/Thigh Support product numbers below, the pads are to be mounted from the Solid Seat Insert.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with. We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Flat Lateral Hip/Thigh Support Pads Mounted From The Seat

Model	Description	Size	Dimension	
			H	D
SHD-NB	Flat Lateral Hip/Thigh Support Pads, pair	X-Small	4	6
SHS-NB	Flat Lateral Hip/Thigh Support Pads, pair	Small	4	9
SHM-NB	Flat Lateral Hip/Thigh Support Pads, pair	Medium	4	12
SHL-NB	Flat Lateral Hip/Thigh Support Pads, pair	Large	4	14
SHO-NB	Flat Lateral Hip/Thigh Support Pads, pair	X-Large	4	16
SHJ-NB	Flat Lateral Hip/Thigh Support Pads, pair	XX-Large	5	16
SHX-NB	Custom Flat Lateral Hip/Thigh Support Pads	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SHL-NB becomes SHL-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


B 1.1.2.4.- Contoured Lateral Hip/Thigh Support Pads Mounted From The Seat

The Contoured Lateral Hip/Thigh Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts as well as keeping the pelvis centered on the seat. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The contouring of the Lateral Hip/Thigh Support Pad will be beneficial by providing a subtle cue to the Trochanteric area to assist in maintaining a properly centered pelvis. The pelvis will be supported and stabilized allowing for more efficient use of the upper and lower extremities. By choosing from the Hip/Thigh Support product numbers below, the pads are to be mounted from the Solid Seat Insert.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Contoured Lateral Hip/Thigh Support Pads Mounted From The Seat

Model	Description	Size	Dimension	
			H	D
SID-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	X-Small	4	6
SIS-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	Small	4	9
SIM-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	Medium	4	12
SIL-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	Large	4	14
SIO-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	X-Large	4	16
SIJ-NB	Contoured Lateral Hip/Thigh Support Pads, contoured at rear, pair	XX-Large	5	16
SIX-NB	Custom Contoured Lateral Hip/Thigh Support Pads, contoured at rear	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SIL-NB becomes SIL-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



B 1.1.2.5.- Flat Lateral Hip/Thigh Support Pads Mounted From Armrests

The Flat Lateral Hip/Thigh Support Pads offered by BioDynamics assist in controlling unwanted lateral weight shifts as well as keeping the pelvis centered on the seat. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity. The pelvic and thigh areas will be supported and stabilized allowing a more efficient use of the upper and lower extremities. By choosing from the Hip/Thigh Support product numbers below, the pads are to be mounted to the armrests.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Flat Lateral Hip/Thigh Support Pads Mounted From Armrests

Model	Description	Size	Dimension	
			H	D
SJD-NB	Flat Lateral Hip/Thigh Support Pads, pair	X-Small	4	8
SJS-NB	Flat Lateral Hip/Thigh Support Pads, pair	Small	4	10
SJM-NB	Flat Lateral Hip/Thigh Support Pads, pair	Medium	4	12
SJL-NB	Flat Lateral Hip/Thigh Support Pads, pair	Large	5	14
SJO-NB	Flat Lateral Hip/Thigh Support Pads, pair	X-Large	5	16
SJJ-NB	Flat Lateral Hip/Thigh Support Pads, pair	XX-Large	5	18
SJX-NB	Custom Flat Lateral Hip/Thigh Support Pads	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.-SJS-NB becomes SJS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed by BioDynamics’ staff, select one of the pre-configured mounting hardware option from section B 1.2.2.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.5.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.1.3.- Lateral Knee Support Pads



B 1.1.3.1.- Flat Lateral Knee Support Pads

The Flat Lateral Knee Support Pads offered by BioDynamics help prevent abduction of the femurs. Properly positioned femurs provide a more stable pelvis allowing for more independent use of the upper extremities. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Flat Lateral Knee Support Pads

Model	Description	Size	Dimension	
			H	D
SLD-NB	Flat Lateral Knee Support Pads, pair	X-Small	3	4
SLS-NB	Flat Lateral Knee Support Pads, pair	Small	4	4
SLM-NB	Flat Lateral Knee Support Pads, pair	Medium	4	5
SLL-NB	Flat Lateral Knee Support Pads, pair	Large	4	6
SLO-NB	Flat Lateral Knee Support Pads, pair	X-Large	5	6
SLX-NB	Custom Flat Lateral Knee Support Pads	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SLS-NB becomes SLS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.1.3.2.- L-Shaped Lateral Knee Support Pads

The L-Shaped Lateral Knee Support Pads offered by BioDynamics help prevent abduction and rotation of the femur. The L-Shaped Pads are helpful when the external rotation exhibited is significant and the additional control of the lower leg will prove helpful. Properly positioned femurs provide a more stable pelvis allowing for more independent use of the upper extremities. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity.

Our pads are made with a Kydex base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Trunk Support Pad he/she may come in contact with.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



L-Shaped Lateral Knee Support Pads

Model	Description	Size	Dimension		
			HR	DT	HF
SND-NB	L-Shaped Lateral Knee Support Pads, pair	X-Small	3	5	5
SNS-NB	L-Shaped Lateral Knee Support Pads, pair	Small	4	5	6
SNM-NB	L-Shaped Lateral Knee Support Pads, pair	Medium	4	6	6
SNL-NB	L-Shaped Lateral Knee Support Pads, pair	Large	4	7	6
SNO-NB	L-Shaped Lateral Knee Support Pads, pair	X-Large	5	7	8
SNX-NB	Custom L-Shaped Lateral Knee Support Pads	Custom	Template or Specs		

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SNS-NB becomes SNS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed in a new or existing BioDynamics Solid Back by our staff, select one of the pre-configured mounting hardware option from section B 1.2.2.1.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.1.3.3.- Round Shaped Lateral Knee Support Pads

The Round Shaped Lateral Knee Support Pads offered by BioDynamics help prevent abduction and rotation of the femur. The Round Shaped Pads are helpful when the external rotation exhibited is significant and the additional control of the lower leg will prove helpful. Properly positioned femurs provide a more stable pelvis allowing for more independent use of the upper extremities. Their use can prevent orthopedic deformities from occurring as well as halt the further deterioration of an existing deformity.

Our pads are made with a H.D.P.E. base covered with a ½” layer of Sunmate foam. The entire perimeter of the pad is also covered with Sunmate foam to offer the consumer protection from all surfaces of the Lateral Knee Support Pad he/she may come in contact with.

This accessory is commonly mounted to the tubing of the footrests or to the side frame of chair. If that is the desire, mounting hardware should be ordered.

We offer a comprehensive selection of standard size pads. If a custom size is needed, select the appropriate product number below and specify size needed.

Pads are sold in pairs.



Round Shaped Lateral Knee Support Pads

Model	Description	Size	Dimension
			Diam.
SKS-NB	Round Shaped Lateral Knee Support Pads, pair	Small	3
SKM-NB	Round Shaped Lateral Knee Support Pads, pair	Medium	4
SKL-NB	Round Shaped Lateral Knee Support Pads, pair	Large	5
SKX-NB	Custom Round Shaped Lateral Knee Support Pads	Custom	Template or Specs

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SKS-NB becomes SKS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed by BioDynamics’ staff, select one of the pre-configured mounting hardware option from section B 1.2.2.3.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.8.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



B 1.1.4.- Medial Knee Support Pads

The line of Medial Knee Support Pads offered by BioDynamics aid in maintaining a neutral position in the lower extremities preventing excessive adduction. The pads are available in 3 standard shapes; L-Shaped, Standard Wedge and Oval.

We begin by using a solid core of wood and surround it with a ½” of Sunmate foam to reduce any pressure areas. The core is t-nutted to match the desired hardware option.

We offer a comprehensive selection of standard size support pads. If a custom size is needed, select the appropriate product number below and specify the dimensions required.

B 1.1.4.1.- Medial Knee Support Pads, Wedge Shaped



Model	Description	Size	Dimension			
			WR	WF	D	H
SMD-NB	Std. Medial Knee Support, wedge shaped, Pad only (pommel)	X-Small	1.5	3	4	3
SMS-NB	Std. Medial Knee Support, wedge shaped, Pad only (pommel)	Small	2	4	4	3
SMM-NB	Std. Medial Knee Support, wedge shaped, Pad only (pommel)	Medium	2	4	5	4
SML-NB	Std. Medial Knee Support, wedge shaped, Pad only (pommel)	Large	3	5	6	4
SMO-NB	Std. Medial Knee Support, wedge shaped, Pad only (pommel)	X-Large	4	6	6	4
SMX-NB	Custom Std. Medial Knee Support, wedge shaped, Pad only (pommel)	Custom	Template or Specs			

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SMS-NB becomes SMS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:



- If the above support pads are being installed by BioDynamics’ staff, select one of the pre-configured mounting hardware option from section B 1.2.2.4.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.7.-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



B 1.1.4.2.- Medial Knee Support Pads, L-Shaped Wedge


Model	Description	Size	Dimension				
			WR	WF	D	HR	HF
SPD-NB	L-Shaped Wedge Medial Knee Support, Pad only (pommel)	X-Small	1.5	3	4	3	5
SPS-NB	L-Shaped Wedge Medial Knee Support, Pad only (pommel)	Small	2	4	4	3	5
SPM-NB	L-Shaped Wedge Medial Knee Support, Pad only (pommel)	Medium	2	4	5	4	6
SPL-NB	L-Shaped Wedge Medial Knee Support, Pad only (pommel)	Large	3	5	6	4	6
SPO-NB	L-Shaped Wedge Medial Knee Support, Pad only (pommel)	X-Large	4	6	6	4	6
SPX-NB	Custom L-Shaped Wedge Medial Knee Support, Pad only (pommel)	Custom	Template or Specs				

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.- SPS-NB becomes SPS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed by BioDynamics' staff, select one of the pre-configured mounting hardware option from section B 1.2.2.4.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.7.-

B 1.1.4.3.- Medial Knee Support Pads, Oval Shaped


Model	Description	Size	Dimension		
			W	D	H
SQD-NB	Oval shaped, Medial Knee Support, Pad only (pommel)	X-Small	2	4	3
SQS-NB	Oval shaped, Medial Knee Support, Pad only (pommel)	Small	3	4	3
SQM-NB	Oval shaped, Medial Knee Support, Pad only (pommel)	Medium	3	5	4
SQL-NB	Oval shaped, Medial Knee Support, Pad only (pommel)	Large	4	6	4
SQX-NB	Custom Oval shaped, Medial Knee Support, Pad only (pommel)	Custom	Template or Specs		

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.- SQS-NB becomes SQS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:



- If the above support pads are being installed by BioDynamics' staff, select one of the pre-configured mounting hardware option from section B 1.2.2.4.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.7.-

“If you don't see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.1.4.4.- Medial Knee Support Pads, L-Shaped Oval


Model	Description	Size	Dimension			
			W	D	HR	HF
SUD-NB	L-Shaped Oval, Medial Knee Support, Pad only (pommel)	X-Small	2	4	3	5
SUS-NB	L-Shaped Oval, Medial Knee Support, Pad only (pommel)	Small	3	4	3	5
SUM-NB	L-Shaped Oval, Medial Knee Support, Pad only (pommel)	Medium	3	5	4	6
SUL-NB	L-Shaped Oval, Medial Knee Support, Pad only (pommel)	Large	4	6	4	6
SUX-NB	Custom L-Shaped Oval, Medial Knee Support, Pad only (pommel)	Custom	Template or Specs			

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- SQS-NB becomes SQS-NE when color Eggplant is desired.

A mounting hardware option needs to be selected:

- If the above support pads are being installed by BioDynamics' staff, select one of the pre-configured mounting hardware option from section B 1.2.2.4.-
- If the above support pads are being installed by RTS, Clinicians, Technician or other select one of the mounting hardware option from section B 1.2.1.7.-


“If you don't see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.2.- Mounting Hardware for Positioning Pads



B 1.2.1.- Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports, Lateral Knee Supports and Medial Knee Supports

The mounting hardware for the Lateral Trunk Supports, Lateral Hip/Thigh Supports, Lateral Knee Supports and the Medial Knee Supports all utilize the most durable materials. BioDynamics' choice of material provides the tensile strength to support significant loads while allowing some flexibility to prevent fractures from occurring in the metal.

For Lateral Trunk, Lateral Hip or Lateral Knee Support Pads, a choice of configuration is necessary to ensure the Pad will extend beyond the cushion sufficiently to allow the pad to provide the clinical benefit necessary. BioDynamics' Medial Knee Support Pads Hardware has significant adjustability to ensure proper cushion clearance. Therefore, only a selection of: Flip Down, Removable or Flip Down/Removable is necessary.


Brackets are sold in pairs.

B 1.2.1.1.- 90 Degree Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

BioDynamics Standard Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports consist of a choice of 90 degree spring steel brackets.

Standard brackets are sold in pairs.

90 Degree Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports






Model	Description	Dimension
 <p>ALA</p>	2" by 1.5" Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Bracket, pair, kit.	(2"x 1.5") by 2.5"wide

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

 <p>ALB</p>	<p>4" by 2" Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Bracket, pair, kit.</p>	<p>(4"x 2") by 2.5"wide</p>
 <p>ALC</p>	<p>4" by 3" Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Bracket, pair, kit.</p>	<p>(4"x 3") by 2.5"wide</p>
 <p>ALD</p>	<p>4" by 4.5" Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Bracket, pair, kit.</p>	<p>(4"x 4.5") by 2.5"wide</p>
 <p>ALE</p>	<p>4" by 6" Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Bracket, pair, kit.</p>	<p>(4"x 6") by 2.5"wide</p>
 <p>ALF</p>	<p>6" by 2" Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Bracket, pair, kit.</p>	<p>(6"x 2") by 2.5"wide</p>

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


B 1.2.1.2.- Offset Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

Offset hardware allows for the Lateral trunk, Lateral Hip/Thigh or Lateral Knee Pads to be positioned within the surface of a solid seat or a solid back. This typically allows for weight gain/loss or seasonal clothing.


When used on either the right or left side of a seat or back, the offset pad often helps to correct a flexible skeletal deformity or accommodate for a fixed skeletal deformity.

A choice of configuration is necessary to allow us to assemble the offset hardware to your specific needs (see drawings). Please note that our offset mounting hardware is fully adjustable and components may be rotated to achieve required dimensions.

Offset brackets are sold in pairs.



AJS



AJL

Offset Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

Model	Description	Offsets
AJS	2ea.(4"x2") + 4ea. (2" x 1.5") Adjustable Offset Bracket Kit for Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Support	Offset Medially: 0.75" to 3" Offset Cushion: 1.5" to 3"
AJL	2ea.(4"x3") + 4ea. (2" x 1.5") Adjustable Offset Bracket Kit for Lateral Trunk, Lateral Hip/Thigh or Lateral Knee Support	Offset Medially: 0.75" to 3" Offset Cushion: 2.75" to 5"

B 1.2.1.3.- Removable Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

The removable hardware options for Lateral Trunk Pads, Lateral Hip/Thigh Pads, and Lateral Knee Pads allow for an easy transfer of the client whether they are a maximum assist or an independent transfer. The pad is removed by pulling a spring loaded button on the bracket. The pad is then easily removed.

The Winter Summer brackets for Lateral Trunk Supports allows for easy lateral adjustment of each pad. With the removal of a set screw on each bracket, the Lateral Trunk Support can be completely removed from the Solid Back Insert.


Removable brackets are sold in pairs.

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



ART to ARM



ARW to ARZ



AOW to AOZ

**Removable Mounting Hardware for Lateral Trunk Supports,
Lateral Hip/Thigh Supports and Lateral Knee Supports**

Model	Description	Size	Dimension
			L
ART	Removable Hardware for Lateral Trunk Supports, kit, pair	Small	3
ARR	Removable Hardware for Lateral Trunk Supports, kit, pair	Medium	5
ARS	Removable Hardware for Lateral Trunk Supports, kit, pair	Large	7
ARH	Removable Hardware for Lateral Hip/Thigh Supports, kit, pair	Small	3
ARI	Removable Hardware for Lateral Hip/Thigh Supports, kit, pair	Medium	5
ARJ	Removable Hardware for Lateral Hip/Thigh Supports, kit, pair	Large	7
ARK	Removable Hardware for Lateral Knee Supports, kit, pair	Small	3
ARL	Removable Hardware for Lateral Knee Supports, kit, pair	Medium	5
ARM	Removable Hardware for Lateral Knee Supports, kit, pair	Large	7
ARW	Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Small	3
ARY	Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Medium	5
ARZ	Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Large	7
AOW	Offset Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Small	3 (offsets 0.75 to 3")
AOY	Offset Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Medium	5 (offsets 0.75 to 3")
AOZ	Offset Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Large	7 (offsets 0.75 to 3")

B 1.2.1.4.- Single Knob and Quick Adjust Assembly for Mounting Hardware of Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

This single knob option offers the consumer/caretaker the ability to slide the pads inwards/outwards and up/down with the turn of a knob. We have found that this option accomplishes the same functionality as the conventional “winter-summer” brackets for a substantially lower cost. This single knob option needs to be ordered with any of our standard brackets that mount Lateral Trunk, Lateral Hip/Thigh, and Lateral Knee Pads.

Knob adjustable mechanisms are sold in pairs.



AKT



AKG

**Single Knob and Quick Adjust Assembly for Mounting Hardware of
Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports**

Model	Description
AKT	Single Knob, Quick Adjust, Assembly to be mounted on a Track Hardware, (winter-summer effect), kit, pair
AKG	Single Knob, Quick Adjust, Assembly to be mounted on a 1/4" groove (1/2" thick), (winter-summer effect), kit, pair

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.2.1.5.- Rubber Coated Metal Straps for Lateral Hip/Thigh Supports and Lateral Knee Supports

BioDynamics' Rubber Coated Straps provide a very quick and easy means to secure Lateral Hip/Thigh and Lateral Knee Pads to the armrest assembly or any other round tubing found on the wheelchair. Pads ordered with this mounting option are shipped with multiple t-nut inserts to ensure a problem free installation.



The Rubber Coated Straps are sold in sets of 4; 2 for each positioning pad.

Rubber Coated Metal Straps for Lateral Hip/Thigh Supports and Lateral Knee Supports

Model	Description
ATA	Rubber Coated Metal Straps for Lateral Hip/Thigh Supports or Lateral Knee Supports, kit of 4ea. for 1" tubing
ATB	Rubber Coated Metal Straps for Lateral Hip/Thigh Supports or Lateral Knee Supports, kit of 4ea. for 7/8" tubing
ATC	Rubber Coated Metal Straps for Lateral Hip/Thigh Supports or Lateral Knee Supports, kit of 4ea. for 3/4" tubing
ATD	Rubber Coated Metal Straps for Lateral Hip/Thigh Supports or Lateral Knee Supports, kit of 4ea. for 5/8" tubing

B 1.2.1.6.- Swing Away Hardware for Lateral Trunk Supports

The Swing Away option is available for Lateral Trunk Supports.. Typically used to effect a lateral transfer, the bracket swings the lateral pad 180 degrees out of position from its functional position to allow for an unencumbered lateral transfer.. Swing Away brackets also allow for the use of a contoured pad when necessary. A contoured pad must be swung away for either a front transfer or a lateral transfer.

The Swing Away feature is easily operated by depressing a spring loaded button on each bracket. When returned to its clinically beneficial position, the bracket will lock the pad in its proper place.

Swing Away brackets are sold in pairs.



AST to ASS



ASO to ASQ



ASW to ASZ



ASI to ASK

Swing Away Hardware for Lateral Trunk Supports

Model	Description	Size	Dimension	
			L	Offsets
AST	Swing Away Hardware for Lateral Trunk Supports, kit, pair	Small	3	n/a
ASR	Swing Away Hardware for Lateral Trunk Supports, kit, pair	Medium	5	n/a
ASS	Swing Away Hardware for Lateral Trunk Supports, kit, pair	Large	7	n/a
ASO	Offset Swing Away Hardware for Lateral Trunk Supports, kit, pair	Small	3	(offsets 0.75 to 3")
ASP	Offset Swing Away Hardware for Lateral Trunk Supports, kit, pair	Medium	5	(offsets 0.75 to 3")
ASQ	Offset Swing Away Hardware for Lateral Trunk Supports, kit, pair	Large	7	(offsets 0.75 to 3")
ASW	Swing Away with Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Small	3	n/a
ASY	Swing Away with Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Medium	5	n/a
ASZ	Swing Away with Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Large	7	n/a
ASI	Offset Swing Away with Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Small	3	(offsets 0.75 to 3")
ASJ	Offset Swing Away with Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Medium	5	(offsets 0.75 to 3")
ASK	Offset Swing Away with Winter Summer Hardware for Lateral Trunk Supports, kit, pair	Large	7	(offsets 0.75 to 3")

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.2.1.7.- Hardware for Medial Knee Supports

All hardware options for mounting a medial knee support pad enable the pad to be moved to allow a smooth transfer of the client. Two main options are offered:


- The Flip Down Bracket for the Medial Knee Support allows the Medial Support Pad to be swung completely under the solid seat to enable a transfer. This Flip Down feature is activated by depressing a spring loaded button on the bracket.
- The removable hardware for the Medial Knee Support allows for the complete removal of the Medial Knee Support to affect a transfer. To remove the pad, the caretaker must unscrew a knob on the underside of the solid seat and the pad will remove. We provide a tether strap to ensure the Medial Knee Pad does not leave the chair during removal.



FMS to FML



RMS to RML



FRS to FRL

Hardware for Medial Knee Supports

Model	Description	Size	Dimension	
			H	D
FMS	Flip Down Hardware for Medial Knee Support	Small	1.5 to 3.5	n/a
FMM	Flip Down Hardware for Medial Knee Support	Medium	2.5 to 5	n/a
FML	Flip Down Hardware for Medial Knee Support	Large	3.5 to 6	n/a
RMS	Removable Hardware for Medial Knee Support	Small	1 to 3.5	0 to 5
RMM	Removable Hardware for Medial Knee Support	Medium	1 to 6	0 to 5
RML	Removable Hardware for Medial Knee Support	Large	1 to 9	0 to 9
FRS	Flip Down and Removable Hardware for Medial Knee Support	Small	1.5 to 3.5	n/a
FRM	Flip Down and Removable Hardware for Medial Knee Support	Medium	2.5 to 5	n/a
FRL	Flip Down and Removable Hardware for Medial Knee Support	Large	3.5 to 6	n/a

B 1.2.1.8.- Hardware for Round Lateral Knee Supports

The Round Shaped Lateral Knee Support Pads are commonly mounted to the tubing of the footrests or to the side frame of chair. We offer a comprehensive selection of standard mounting hardware for these pads.

Hardware for Round Lateral Knee Supports

Model	Description	
 # ATM & ATN	ATM	Rubber Coated Metal Straps for Round Lateral Knee Supports, kit of 2ea. for 1" tubing
	ATN	Rubber Coated Metal Straps for Round Lateral Knee Supports, kit of 2ea. for 7/8" tubing

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

 <p># ATO & ATP</p>	ATO	Rubber Coated Metal Straps for Round Lateral Knee Supports, kit of 2ea. for 1" tubing and 2ea. Bracket HJJ
	ATP	Rubber Coated Metal Straps for Round Lateral Knee Supports, kit of 2ea. for 7/8" tubing and 2ea. Bracket HJJ
 <p># ATQ & ATR</p>	ATQ	2ea. (1/4"x20) 1.75" Button Head Bolt (mounts to the pivot point of ELR's)
	ATR	2ea. (5/16"x20) 1.75" Button Head Bolt (mounts to the pivot point of ELR's)

B 1.2.2.- Pre-Configured Mounting Hardware for Positioning Pads

B 1.2.2.1.- Pre-Configured Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

When positioning pads are chosen as part of a seating system assembled at BioDynamics, maximum growth and adjustability are obtained by selecting a pre-configured mounting option from the table below.

The mounting hardware for the Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports all utilize the most durable materials. BioDynamics' choice of material provides the tensile strength to support significant loads while allowing some flexibility to prevent fractures from occurring in the metal.

Brackets are sold in pairs.

Pre-Configured Mounting Hardware for Lateral Trunk Supports, Lateral Hip/Thigh Supports and Lateral Knee Supports

Model	Description
 <p>LBB</p>	Straight Bracket (90 degree Spring Steel).

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

 <p style="text-align: center;">LBK</p>	<p>Straight Bracket w/ Single Knob.</p>
 <p style="text-align: center;">LSB</p>	<p>Swing Away Bracket.</p>
 <p style="text-align: center;">LSW</p>	<p>Swing Away w/ Winter Summer Bracket.</p>
 <p style="text-align: center;">LSK</p>	<p>Swing Away w/ Single Knob Bracket.</p>
 <p style="text-align: center;">LWB</p>	<p>Winter Summer w/Straight Bracket (90 degree).</p>

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

 <p style="text-align: center;">LOF</p>	<p>Removable Bracket for Lateral Trunk, Lateral Hip/Thigh and Lateral Knee Support.</p>
 <p style="text-align: center;">LOB</p>	<p>Offset Bracket (Adjusts from .75" to 3" on each side).</p>
 <p style="text-align: center;">LOK</p>	<p>Offset Bracket w/ Single Knob.</p>
 <p style="text-align: center;">LOS</p>	<p>Offset Swing Away Bracket.</p>
 <p style="text-align: center;">LOR</p>	<p>Offset Swing Away w/ Winter Summer Bracket.</p>

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

 LOT	Offset Swing Away w/ Single Knob Bracket.
 LOW	Offset Winter Summer.
LOX	Custom Mounting Hardware

B 1.2.2.2.- Pre-Configured Mounting Hardware for Lateral Hip/Thigh Supports secured from the Armrest

When Lateral Hip/Thigh Supports pads are chosen to be secured from the armrest and mounted by BioDynamics, select an option from the table below.

Brackets are sold in pairs.



#ATH

Pre-Configured Mounting Hardware for Lateral Hip/Thigh Supports secured from the Armrest

Model	Description
ATH	Mounting Hardware for Hip/Thigh to Tubular Armrest Frame by Rubber Coated Metal Straps.
ATX	Custom Mounting Hardware for Hip/Thigh to Non-tubular Armrest Frame

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 1.2.2.3.- Pre-Configured Mounting Hardware for Round Lateral Knee Supports

When Round Lateral Knee Supports pads are chosen to be mounted by BioDynamics, select an option from the table below.

Brackets are sold in pairs.



#AUL



#AUE



#AUF

Pre-Configured Mounting Hardware for Round Lateral Knee Supports

Model	Description
AUL	Mounting Hardware for Positioning Pads to Footrest
AUE	Mounting Hardware for Positioning Pads to ELR
AUF	Mounting Hardware for Positioning Pads to Side Frame
AUS	Mounting Hardware for Positioning Pads to front of Seat Pan
AUX	Custom Mounting Hardware for Positioning Pads

B 1.2.2.4.- Pre-Configured Mounting Hardware for Medial Knee Supports

When Medial Knee Supports pads are chosen to be mounted by BioDynamics, select an option from the table below.



#AWF



#AWR



#AWB

Pre-Configured Mounting Hardware for Medial Knee Supports

Model	Description
AWF	Flip Down Hardware for Medial Knee Support
AWR	Removable Hardware for Medial Knee Support
AWB	Flip Down and Removable Hardware for Medial Knee Support
AWX	Custom Mounting of Medial Knee Support

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.- Positioning Belts, Straps and Pads

Positioning Belts, Straps and Pads found in this section enhance the positioning of the upper trunk and lower extremities. Many of the Positioning Belts and Straps are available with pads. Those without can be retrofitted with our padded sleeve options. All choices come with the appropriate mounting hardware.

B 2.1- Anterior Trunk Supports



The BioForm Line of Anterior Trunk Supports are available in Neoprene (Rubatex), Reverse Dartex, Vinyl (Naugahyde), and Ballistic Nylon. Neoprene is effective when control is needed but a stretch component would prove beneficial. Reverse Dartex will offer a more rigid support than Neoprene while still having some elasticity, as well as providing an easy to clean surface. Vinyl (Naugahyde) will offer the most rigid support as well as a surface which handles fluids and cleanings most easily. Ballistic Nylon is the preferred material when the consumer will try to orally or otherwise try to destroy the harness.

All styles below are available in standard sizes from XS through XL. A customization to the size and/or shape is available by selecting the appropriate model #.

“If you don’t see it... Contact Us”


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

- **BioForm Standard Anterior Trunk Supports** are an effective way to help maintain a properly positioned trunk when significant assistance is needed. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm Anterior Trunk Supports with Snap Closures** assist in easier/quicker transfers from the wheelchair. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm “H” type Anterior Trunk Supports** are effective when less than maximum assistance is needed to maintain upper trunk position and quick and easy transfer is a primary goal. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm Low Profile Anterior Trunk Supports** are effective when significant assistance is needed to properly maintain upper trunk position and a g-tube needs to be accommodated for. This product may also prove helpful when used with the adolescent to fully matured female shape. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm 6 Point Trunk Support** is effective for the client who needs **maximum** trunk support: Anterior Trunk, pelvis and shoulders. Comes complete with 6 securing straps, end fittings and all mounting hardware.

B 2.1.1.- BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support



Model	Description	Size	Dimension					
			A	B	C	D	E	F
NRT-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
NRD-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
NRS-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	Small	17	9	10	5	2	5
NRM-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
NRL-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
NRO-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
NRX-NB	BioForm Stretch (Black Neoprene) Custom Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NRD-NB becomes NRD-NE when color Eggplant is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.2.- BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
NLD-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	X-Small	14	8	8.5	4	2	1
NLS-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	Small	17	9	10	5	2	1
NLM-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	Medium	20	10	11.5	6	2	1.25
NLL-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	Large	22	11	12	7	2.5	1.5
NLX-NB	BioForm Stretch (Black Neoprene) Custom Low-Profile Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NLD-NB becomes NLD-NE when color Eggplant is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.3.- BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure


Model	Description	Size	Dimension					
			A	B	C	D	E	F
NST-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	XX-Small	10	5	5	4	1.75	1.25
NSD-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	X-Small	14	8	8.5	4	2	4
NSS-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	Small	17	9	10	5	2	5
NSM-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	Medium	20	10	11.5	6	2	6
NSL-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	Large	22	11	12	7	2.5	7
NSO-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	X-Large	24	12	14	8	2.5	8
NSX-NB	BioForm Stretch (Black Neoprene) Custom Anterior Trunk Support, Snap Closure	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NSD-NB becomes NSD-NE when color Eggplant is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.4.- BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support


Model	Description	Size	Dimension				
			A	B	C	D	E
NHS-NB	BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support	Small	10	14	6	8	2
NHM-NB	BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support	Medium	13	18	8	10	2.25
NHL-NB	BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support	Large	16	21	10	12	2.5
NHX-NB	BioForm Stretch (Black Neoprene) Custom “H”-Style Anterior Trunk Support	Custom	Template or Specs				

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NHD-NB becomes NHD-NE when color Eggplant is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.5.- BioForm Stretch (Black Neoprene) 6 Point Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
NMT-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	XX-Small	10	5	5	4	1.75	1.25
NMD-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	X-Small	14	8	8.5	4	2	4
NMS-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	Small	17	9	10	5	2	5
NMM-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	Medium	20	10	11.5	6	2	6
NML-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	Large	22	11	12	7	2.5	7
NMO-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	X-Large	24	12	14	8	2.5	8
NMX-NB	BioForm Stretch (Black Neoprene) Custom 6 Point Trunk Support	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NMD-NB becomes NMD-NE when color Eggplant is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.6.- BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
DRT	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
DRD	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
DRS	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	Small	17	9	10	5	2	5
DRM	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
DRL	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
DRO	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
DRX	BioForm Semi-Stretch (Black Reverse Dartex) Custom Anterior Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.7.- BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support



Model	Description	Size	Dimension					
			A	B	C	D	E	F
DLD	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	X-Small	14	8	8.5	4	2	1
DLS	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	Small	17	9	10	5	2	1
DLM	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	Medium	20	10	11.5	6	2	1.25
DLL	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	Large	22	11	12	7	2.5	1.5
DLX	BioForm Semi-Stretch (Black Reverse Dartex) Custom Low-Profile Anterior Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

B 2.1.8.- BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure


Model	Description	Size	Dimension					
			A	B	C	D	E	F
DST	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	XX-Small	10	5	5	4	1.75	1.25
DSD	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	X-Small	14	8	8.5	4	2	4
DSS	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	Small	17	9	10	5	2	5
DSM	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	Medium	20	10	11.5	6	2	6
DSL	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	Large	22	11	12	7	2.5	7
DSO	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	X-Large	24	12	14	8	2.5	8
DSX	BioForm Semi-Stretch (Black Reverse Dartex) Custom Anterior Trunk Support, Snap Closure	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.9.- BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support


Model	Description	Size	Dimension				
			A	B	C	D	E
DHS	BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support	Small	10	14	6	8	2
DHM	BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support	Medium	13	18	8	10	2.25
DHL	BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support	Large	16	21	10	12	2.5
DHX	BioForm Semi-Stretch (Black Reverse Dartex) Custom “H”-Style Anterior Trunk Support	Custom	Template or Specs				


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.10.- BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
DMT	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	XX-Small	10	5	5	4	1.75	1.25
DMD	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	X-Small	14	8	8.5	4	2	4
DMS	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	Small	17	9	10	5	2	5
DMM	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	Medium	20	10	11.5	6	2	6
DML	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	Large	22	11	12	7	2.5	7
DMO	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	X-Large	24	12	14	8	2.5	8
DMX	BioForm Semi-Stretch (Black Reverse Dartex) Custom 6 Point Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.11.- BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VRT-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
VRD-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
VRS-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	Small	17	9	10	5	2	5
VRM-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
VRL-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
VRO-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
VRX-V393	BioForm Non-Stretch (Black Vinyl) Custom Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VRD-V393 becomes VRD-V322 when color Sun Yellow is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.12.- BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VLD-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	X-Small	14	8	8.5	4	2	1
VLS-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	Small	17	9	10	5	2	1
VLM-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	Medium	20	10	11.5	6	2	1.25
VLL-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	Large	22	11	12	7	2.5	1.5
VLX-V393	BioForm Non-Stretch (Black Vinyl) Custom Low-Profile Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VLD-V393 becomes VLD-V322 when color Sun Yellow is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.13.- BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VST-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	XX-Small	10	5	5	4	1.75	1.25
VSD-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	X-Small	14	8	8.5	4	2	4
VSS-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	Small	17	9	10	5	2	5
VSM-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	Medium	20	10	11.5	6	2	6
VSL-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	Large	22	11	12	7	2.5	7
VSO-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	X-Large	24	12	14	8	2.5	8
VSX-V393	BioForm Non-Stretch (Black Vinyl) Custom Anterior Trunk Support, Snap Closure	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VSD-V393 becomes VSD-V322 when color Sun Yellow is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.14.- BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support


Model	Description	Size	Dimension				
			A	B	C	D	E
VHS-V393	BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support	Small	10	14	6	8	2
VHM-V393	BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support	Medium	13	18	8	10	2.25
VHL-V393	BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support	Large	16	21	10	12	2.5
VHX-V393	BioForm Non-Stretch (Black Vinyl) Custom “H”-Style Anterior Trunk Support	Custom	Template or Specs				

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VHD-V393 becomes VHD-V322 when color Sun Yellow is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.15.- BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VMT-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	XX-Small	10	5	5	4	1.75	1.25
VMD-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	X-Small	14	8	8.5	4	2	4
VMS-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	Small	17	9	10	5	2	5
VMM-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	Medium	20	10	11.5	6	2	6
VML-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	Large	22	11	12	7	2.5	7
VMO-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	X-Large	24	12	14	8	2.5	8
VMX-V393	BioForm Non-Stretch (Black Vinyl) Custom 6 Point Trunk Support	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VMD-V393 becomes VMD-V322 when color Sun Yellow is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-


“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.16.- BioForm Black Ballistic Nylon Standard Anterior Trunk Support




Model	Description	Size	Dimension					
			A	B	C	D	E	F
VNT	BioForm Ballistic Nylon Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
VND	BioForm Ballistic Nylon Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
VNS	BioForm Ballistic Nylon Standard Anterior Trunk Support	Small	17	9	10	5	2	5
VNM	BioForm Ballistic Nylon Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
VNL	BioForm Ballistic Nylon Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
VNO	BioForm Ballistic Nylon Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
VNX	BioForm Ballistic Nylon Custom Anterior Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.2.- BioForm Abdominal Binders

The BioForm Abdominal Binders provide a crafty solution for proper positioning of a clients' mid section. It is constructed from a ½” piece of Sunmate foam covered in Neoprene with four straps attached to the body of the Abdominal Binder. Straps are equipped with Tri-Glides (for easy adjustment) and End fittings for easy mounting. Available in 3 sizes.

This option is a sensible choice when presented with the problem of a lordotic midsection that you would like to hold rearwards in the wheelchair.



BioForm Abdominal Binders

Model	Description	Size	Dimension	
			W	H
NAS-NB	BioForm Abdominal Binder, Four 1” Wide Straps, Fastex Buckle Closure, End fittings	Small	6	4
NAM-NB	BioForm Abdominal Binder, Four 1” Wide Straps, Fastex Buckle Closure, End fittings	Medium	9	5.5
NAL-NB	BioForm Abdominal Binder, Four 1” Wide Straps, Fastex Buckle Closure, End fittings	Large	12	7
NAX-NB	BioForm Custom Abdominal Binder	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- NAS-NB becomes NAS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.3.- BioForm Positioning Belts

The BioForm positioning belts present a perfect complement to our seating systems. BioForm Positioning Belts are typically used at a 45 degree angle to properly position the pelvis or at a 90 degree angle across the femurs to maintain the leg contact with the seating surface.

The belts are available in 1", 1½" and 2" webbing widths and come with a choice of 3 closures; d-ring/velcro, fastex buckle and auto style buckle. Every belt comes standard with a tri-glide (for easy adjustment) and a d-ring /end fitting unit for easy mounting.

A padded belt can be created by choosing one of our padded sleeve options.

B 2.3.1.- BioForm Positioning Belts, Velcro Closure

	Model	Description	Dimension	
			W	L
	PBB	BioForm Positioning Belt, 1" webbing, D-Ring and Velcro Closure, End fittings	1	50
	PBC	BioForm Positioning Belt, 1 1/2" webbing, D-Ring and Velcro Closure, End fittings	1.5	60
	PBD	BioForm Positioning Belt, 2" webbing, D-Ring and Velcro Closure, End fittings	2	70
PBX	BioForm Custom Positioning Belt, webbing, D-Ring and Velcro Closure, End fittings	Template or Specs		

B 2.3.2.- BioForm Positioning Belts, Fastex Buckle Closure

	Model	Description	Dimension	
			W	L
	PEB	BioForm Positioning Belt, 1" webbing, Fastex Buckle Closure, End fittings	1	50
	PEC	BioForm Positioning Belt, 1 1/2" webbing, Fastex Buckle Closure, End fittings	1.5	60
	PED	BioForm Positioning Belt, 2" webbing, Fastex Buckle Closure, End fittings	2	70
PEX	BioForm Custom Positioning Belt, webbing, Fastex Buckle Closure, End fittings	Template or Specs		

B 2.3.3.- BioForm Positioning Belts, Auto Style Buckle Closure

	Model	Description	Dimension	
			W	L
	PFB	BioForm Positioning Belt, 1" webbing, Auto Style Buckle Closure, End fittings	1	50
	PFC	BioForm Positioning Belt, 1 1/2" webbing, Auto Style Buckle Closure, End fittings	1.5	60
	PFD	BioForm Positioning Belt, 2" webbing, Auto Style Buckle Closure, End fittings	2	70
PFX	BioForm Custom Positioning Belt, webbing, Auto Style Buckle Closure, End fittings	Template or Specs		

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.4.- BioForm Positioning Belts with Split Pads

The BioForm Positioning Belts are also offered with “split pads” secured to both sides of the positioning belt. This option is effective for accurate protection of the A.S.I.S. bones as the pads can not slide off of the bony prominences.

The belts are available in 1”, 1½” and 2” webbing widths and come with a choice of 3 closures; d-ring/velcro, fastex buckle and auto style buckle. Every belt comes standard with a tri-glide (for easy adjustment) and a d-ring /end fitting unit for easy mounting.

B 2.4.1.- BioForm Positioning Belts, Split Pads, 2 each, Fastex Buckle Closure

	Model	Description	Dimension			
			Belt W	Belt L	Pad W	Pad L
	PIS	BioForm Positioning Belt, 1” wide webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	1	50	1.75	4
	PIT	BioForm Positioning Belt, 1½” wide webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	1.5	60	2.5	6
	PIU	BioForm Positioning Belt, 2” wide webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	2	70	4	8
	PIX	BioForm Custom Positioning Belt, webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	Template or Specs			

B 2.4.2.- BioForm Positioning Belts, Split Pads, 2 each, Auto-Style Buckle Closure

	Model	Description	Dimension			
			Belt W	Belt L	Pad W	Pad L
	PJV	BioForm Positioning Belt, 1” wide webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1	50	1.75	4
	PJW	BioForm Positioning Belt, 1½” wide webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1.5	60	2.5	6
	PJY	BioForm Positioning Belt, 2” wide webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	2	70	4	8
	PJX	BioForm Custom Positioning Belt, webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	Template or Specs			

B 2.4.3.- BioForm 4-Points Positioning Belts, Split Pads, 2 each, Auto-Style Buckle Closure

	Model	Description	Dimension			
			Belt W	Belt L	Pad W	Pad L
	PQV	BioForm 4-Points Positioning Belt, 1” wide webbing (¾” for anchoring), Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1	50	1.75	4
	PQW	BioForm 4-Points Positioning Belt, 1½” wide webbing (1” for anchoring), Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1.5	60	2.5	6
	PQY	BioForm 4-Points Positioning Belt, 2” wide webbing (1” for anchoring), Split Pads (2ea), Auto-Style Buckle Closure, End fittings	2	70	4	8
	PQX	BioForm Custom 4-Points Positioning Belt, webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	Template or Specs			

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.5.- BioForm Chest Positioning Belts

The BioForm Chest Positioning Belts are beneficial for the consumer who needs assistance in keeping their upper trunk in contact to the solid back. These Positioning Belts can be installed to the Lateral Trunk Supports, to the Solid Back Inserts or to the wheelchair frame.

B 2.5.1.- BioForm Non-Stretch Chest Positioning Belts

The BioForm Non-Stretch Chest Positioning Belts are available in 1½” and 2” webbing widths. Both come with d-ring and velcro closure to allow for greater leverage when tightening. One side of the strap is placed in front of the client and passed through the d-ring on the other strap and brought back upon itself and secured via velcro. Positioning Belts have tri-glides for easy adjustment and end fittings for convenient mounting.

A padded Chest Positioning Belt can be created by adding one of the padded sleeve options.

BioForm Non-Stretch Chest Positioning Belts

	Model	Description	Dimension	
			W	L
	PKB	BioForm Chest Positioning Belt, 1½” wide webbing, D-Ring and Velcro Closure, End fittings	1.5	60
	PKC	BioForm Chest Positioning Belt, 2” wide webbing, D-Ring and Velcro Closure, End fittings	2	70
	PKX	BioForm Custom Chest Positioning Belt, webbing, D-Ring and Velcro Closure, End fittings	Template or Specs	

B 2.5.2.- BioForm Stretch (Black Neoprene) Chest Positioning Belts

The BioForm Stretch (Black Neoprene) Chest Positioning Belts are available with a pad of 8”, 11”, 14”, 17” and 20” in length. They come with velcro closure to allow for a quick setup. Each strap, short piece of webbing sewn to the Neoprene pads, overlap each other and attach via velcro in front of the client. Positioning Belts have tri-glides for easy adjustment and end fittings for convenient mounting. Available in x-small, small, medium, large and x-large. Custom sizes available by selecting the appropriate model #.

BioForm Stretch Chest Positioning Belts

	Model	Description	Size	Dimension (overlap)	
				W	L
	PDD-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 1.5” End fittings	X-Small	2.5	6
	PDS-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	Small	4	9
	PDM-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	Medium	4	12
	PDL-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	Large	4	15
	PDO-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	X-Large	5	20
	PDX-NB	BioForm Custom Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, End fittings	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- PDD-NB becomes PDD-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.6.- BioForm Calf Straps

The BioForm Calf Strap is a simple solution to assist in holding the lower extremities from slipping rearwards between the legrests. The strap is 2” wide webbing with a velcro closure. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.

The Calf Strap can be made padded by adding one of the BioForm Padded Sleeve options.

BioForm Calf Straps

	Model	Description	Size	Dimension	
				W	L
	PCS	BioForm Calf Strap, 2” wide webbing with a Velcro closure	Small		30
	PCM	BioForm Calf Strap, 2” wide webbing with a Velcro closure	Medium		35
	PCL	BioForm Calf Strap, 2” wide webbing with a Velcro closure	Large		40
	PCX	BioForm Custom Calf Strap, webbing with a Velcro closure	Custom	Template or Specs	

B 2.7.- BioForm Thigh Straps

The BioForm Thigh Straps are typically chosen when a belt positioned at 90 degrees to the thighs is useful in keeping the femurs in contact with the seat cushion. The straps are available in 1½” and 2” wide webbing with D-Ring and Velcro Closure.

The Thigh Strap can be made padded by adding one of the BioForm Padded Sleeve options.

BioForm Thigh Straps

	Model	Description	Dimension	
			W	L
	PTB	BioForm Thigh Strap, 1½” wide webbing, D-Ring and Velcro Closure, End fittings	1.5	60
	PTC	BioForm Thigh Strap, 2” wide webbing, D-Ring and Velcro Closure, End fittings	2	70
	PTX	BioForm Custom Thigh Strap, webbing, D-Ring and Velcro Closure, End fittings	Template or Specs	

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.8.- BioForm Padded Sleeves



The BioForm Padded Sleeves are the perfect selection when any belt currently in use produces excessive pressure areas. They can be retrofitted to any positioning belt.

The BioForm line of Padded Sleeves are made from ½” of Sunmate foam covered in Neoprene (Rubatex). Two attaching loops are sewn onto pad. One side is fixed to allow belt or strap to pass through and the other is a Velcro overlap to allow for quick and easy removal so belt or strap can be easily undone. This securing method assures caregivers that the pad will stay attached to belt or strap while still easily removable for transfers. Pads are available to fit our 1 ½” belts and our 2” belts and come in 10”, 12” and 14” lengths. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.

B 2.8.1.- BioForm Black Padded Sleeves For 2” Belts

Model	Description	Size	Dimension	
			L	W
PGS-NB	BioForm Black Padded Sleeve for 2” Belt	Small	10”	
PGM-NB	BioForm Black Padded Sleeve for 2” Belt	Medium	12.5”	
PGL-NB	BioForm Black Padded Sleeve for 2” Belt	Large	15”	
PGX-NB	BioForm Custom Black Padded Sleeve for 2” Belt	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- PGM-NB becomes PGM-NE when color Eggplant is desired.

B 2.8.2.- BioForm Black Padded Sleeves For 1½” Belts

Model	Description	Size	Dimension	
			L	W
PHS-NB	BioForm Black Padded Sleeve for 1½” Belt	Small	10”	
PHM-NB	BioForm Black Padded Sleeve for 1½” Belt	Medium	12.5”	
PHL-NB	BioForm Black Padded Sleeve for 1½” Belt	Large	15”	
PHX-NB	BioForm Custom Black Padded Sleeve for 1½” Belt	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- PHM-NB becomes PHM-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.- Lower Extremity Supports

The Lower Extremity Supports offered in this section enable proper positioning of the lower extremities. Every selection comes with a mounting hardware option.

B 3.1.- BioForm Foot Positioners



The BioForm line of foot and ankle positioners include: Ankle Positioners, Footplate Extenders, Ankle Positioners with adjustable Plate and Foot Positioners,

- **BioForm Ankle Positioners** are a 3 point ankle cuff forming a cup around the ankle to maintain proper positioning of the foot/ankle on a footplate. They can be mounted to any slotted footplate or our BioForm Multi-Slots Footplate Extenders. All mounting hardware is included. Fastex Buckle closure. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.
- **BioForm Multi-Slots Footplate Extenders** allow for the foot to be positioned in a location off of the footplate; anterior, posterior or angled. To these multi-slots extenders, our foot and ankle positioning accessories can be mounted in the multiple mounting slots available. All mounting hardware is included. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.
- **BioForm Ankle Positioners with adjustable Plate** includes a BioForm Footplate Extender and a BioForm Ankle Positioner allowing for a proper positioning of the foot/ankle in a location of the footplate. All mounting hardware is included. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.
- **BioForm Foot Positioners** provide for maximum positioning capability of the foot and ankle. They consist of a BioForm Multi-Slots Footplate Extender, BioForm Ankle Positioner and a BioForm Padded Ankle/Toe Strap. All mounting hardware is included. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model.

B 3.1.1.- BioForm Ankle Positioners

	Model	Description	Size	Fits Ankle Circumference	
				From	To
	FAS-NB	BioForm Ankle Positioner, pair	Small	5	9
	FAM-NB	BioForm Ankle Positioner, pair	Medium	7.5	12
	FAL-NB	BioForm Ankle Positioner, pair	Large	10	15
FAX-NB	BioForm Custom Ankle Positioner	Custom	Template or Specs		


If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- FAS-NB becomes FAS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.1.2.- BioForm Multi-Slot Footplate Extenders


Model	Description	Size	Dimension	
			W	L
FXS	BioForm Multi-Slot Footplate Extenders, pair	Small	4	8
FXM	BioForm Multi-Slot Footplate Extenders, pair	Medium	5	9.5
FXL	BioForm Multi-Slot Footplate Extenders, pair	Large	5	11
FXX	BioForm Custom Multi-Slot Footplate Extenders	Custom	Template or Specs	

B 3.1.3.- BioForm Ankle Positioner with Adjustable Plate


Model	Description	Size	Fits Ankle Circumference		Dimension (Plate)	
			From	To	W	L
FBS-NB	BioForm Ankle Positioner With Adjustable Plate, pair	Small	5	9	4	8
FBM-NB	BioForm Ankle Positioner With Adjustable Plate, pair	Medium	7.5	12	5	9.5
FBL-NB	BioForm Ankle Positioner With Adjustable Plate, pair	Large	10	15	5	11
FBX-NB	BioForm Custom Ankle Positioner With Adjustable Plate	Custom	Template or Specs			

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- FBS-NB becomes FBS-NE when color Eggplant is desired.

B 3.1.4.- BioForm Foot Positioner


Model	Description	Size	Fits Ankle Circumference		Dimension (Plate)	
			From	To	W	L
FPS-NB	BioForm Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender) , pair	Small	5	9	4	8
FPM-NB	BioForm Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender) , pair	Medium	7.5	12	5	9.5
FPL-NB	BioForm Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender) , pair	Large	10	15	5	11
FPX-NB	BioForm Custom Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender)	Custom	Template or Specs			

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- FPS-NB becomes FPS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.2.- BioForm Footplates


B 3.2.1.- BioForm One Piece, Footplates

BioDynamics manufactures One Piece, Footplates that are shipped with pre-drilled holes (countersunk) and mounting hardware for easy mounting in the appropriate position. Our One Piece, Footplates are made with a 1/2" High Density Polyethylene (H.D.P.E.) base. The padded version has a 1/4" of H.D.P.E., summate foam and Vinyl (Naugahyde) covered module secured to the base with hardware.

Padded and non-padded One Piece, Footplates are available with fixed hardware or flip up hardware.. The flip up hardware allows the footrests to be swung away/removed when necessary. When the footrests are put back on the chair, the One Piece, Footplate swings down to the opposite footplate from where it is attached and has an indexing pin that will set into a hole the RTS will drill into that footplate. This makes the One Piece, Footplate more secure when in operation.


Our standard sizes are available to fit 14", 16", 18" and 20" wheelchairs. If a custom size is needed, select the appropriate model#.

B 3.2.1.1.- BioForm One Piece, Fixed Footplate



Model	Description	Dimension	
		W	D
CFS	BioForm One Piece, Fixed Footplate, for 14" Wheelchairs	11.5	8
CFM	BioForm One Piece, Fixed Footplate, for 16" Wheelchairs	13.5	9
CFL	BioForm One Piece, Fixed Footplate, for 18" Wheelchairs	15.5	10
CFO	BioForm One Piece, Fixed Footplate, for 20" Wheelchairs	17.5	10
CFX	BioForm Custom One Piece, Fixed Footplate	Template or Specs	

B 3.2.1.2.- BioForm Padded One Piece, Fixed Footplate



Model	Description	Dimension	
		W	D
CGS-V393	BioForm Padded One Piece, Fixed Footplate, for 14" Wheelchairs	11.5	8
CGM-V393	BioForm Padded One Piece, Fixed Footplate, for 16" Wheelchairs	13.5	9
CGL-V393	BioForm Padded One Piece, Fixed Footplate, for 18" Wheelchairs	15.5	10
CGO-V393	BioForm Padded One Piece, Fixed Footplate, for 20" Wheelchairs	17.5	10
CGX-V393	BioForm Padded Custom One Piece, Fixed Footplate	Template or Specs	


If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- CGS-V393 becomes CGS-V354 when Vinyl Royal Blue is desired.

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.2.1.3.- BioForm One Piece, Flip-Up Footplate


Model	Description	Dimension	
		W	D
CES	BioForm One Piece, Flip-Up Footplate, Self Centering Locking, for 14" Wheelchairs	11.5	8
CEM	BioForm One Piece, Flip-Up Footplate, Self Centering Locking,, for 16" Wheelchairs	13.5	9
CEL	BioForm One Piece, Flip-Up Footplate, Self Centering Locking, for 18" Wheelchairs	15.5	10
CEO	BioForm One Piece, Flip-Up Footplate, Self Centering Locking, for 20" Wheelchairs	17.5	10
CEX	BioForm Custom One Piece, Flip-Up Footplate, Self Centering Locking,	Template or Specs	

B 3.2.1.4.- BioForm Padded One Piece, Flip-Up Footplate


Model	Description	Dimension	
		W	D
CDS-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 14" Wheelchairs	11.5	8
CDM-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 16" Wheelchairs	13.5	9
CDL-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 18" Wheelchairs	15.5	10
CDO-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 20" Wheelchairs	17.5	10
CDX-V393	BioForm Custom Padded One Piece, Flip-Up Footplate, Self Centering Locking,	Template or Specs	

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- CDS-V393 becomes CDS-V354 when Vinyl Royal Blue is desired.

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.2.2.- BioForm Footplates Extenders

BioDynamics manufactures Footplates Extenders that are shipped with pre-drilled holes (countersunk) and mounting hardware for easy mounting in the appropriate position. Our Footplates Extenders are made with a 1/2" High Density Polyethylene (H.D.P.E.) base. The padded version has a 1/4" of H.D.P.E., summate foam and Vinyl (Naugahyde) covered module secured to the base with hardware.



Our standard sizes are available to fit 14", 16", 18" and 20" wheelchairs. If a custom size is needed, select the appropriate model#.

B 3.2.2.1.- BioForm Footplates Extender



Model	Description	Dimension	
		W	D
CIS	BioForm Footplates Extender, for 14" Wheelchairs, pair.	5.5	9
CIM	BioForm Footplates Extender, for 16" Wheelchairs, pair.	6.5	10
CIL	BioForm Footplates Extender, for 18" Wheelchairs, pair.	7.5	11
CIO	BioForm Footplates Extender, for 20" Wheelchairs, pair.	8.5	12
CIX	BioForm Custom Footplates Extender	Template or Specs	

B 3.2.2.2.- BioForm Padded Footplates Extender



Model	Description	Dimension	
		W	D
CJS-V393	BioForm Padded Footplates Extender, for 14" Wheelchairs, pair.	5.5	9
CJM-V393	BioForm Padded Footplates Extender, for 16" Wheelchairs, pair.	6.5	10
CJL-V393	BioForm Padded Footplates Extender, for 18" Wheelchairs, pair.	7.5	11
CJO-V393	BioForm Padded Footplates Extender, for 20" Wheelchairs, pair.	8.5	12
CJX-V393	BioForm Custom Padded Footplates Extender	Template or Specs	

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- CJS-V393 becomes CJS-V354 when Vinyl Royal Blue is desired.

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.2.3.- BioForm Footplate Pads

B 3.2.3.1.- BioForm Soft Footplate Pads

BioDynamics manufactures BioForm Soft Footplate Pads that are designed to protect the consumers feet from the footplates. The Footplate Pads have a pocket sewn onto the underside of the pad which enables the entire pad to slip onto the footplate. The pad is then further secured to the footplate with a velcro strap, also, sewn onto the pad. The pads are made of ½” of Sunmate foam covered in Neoprene.

Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#. Sold in pairs.



BioForm Soft Footplate Pads

Model	Description	Dimension	
		W	D
CMS-NB	BioForm Soft Footplate Pads, for 14” Wheelchairs, pair	11.5	8
CMM-NB	BioForm Soft Footplate Pads, for 16” Wheelchairs, pair	13.5	9
CML-NB	BioForm Soft Footplate Pads, for 18” Wheelchairs, pair	15.5	10
CMO-NB	BioForm Soft Footplate Pads, for 20” Wheelchairs, pair	17.5	10
CMX-NB	BioForm Custom Soft Footplate Pads	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- CMS-NB becomes CMS-NE when color Eggplant is desired.

B 3.2.3.2.- BioForm Footplate Pads

BioDynamics manufactures Footplate Pads that are designed to protect the consumers feet from the footplates. These Footplate Pads are secured to the footplate by inserting a pair of bolts through the underside of the footplate and into the T-nuts embedded into the plastic that is on the underside of the Footplate Pad. The pads are made of ½” of Sunmate foam covered in Neoprene which is secured to a ¼” plastic (HDPE) base.

Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#. Sold in pairs.



BioForm Footplate Pads

Model	Description	Size	Dimension	
			W	D
CNS-NB	BioForm Footplate Pads, pair	Small	4	6.5
CNM-NB	BioForm Footplate Pads, pair	Medium	5	6.5
CNL-NB	BioForm Footplate Pads, pair	Large	6	8
CNO-NB	BioForm Footplate Pads, pair	X-Large	7	8
CNX-NB	BioForm Custom Footplate Pads	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- CNS-NB becomes CNS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.3.- BioForm Sandals and Straps



Our BioForm Sandals securely hold the feet in proper position utilizing either: an adjustable Fastex Buckle strap, adjustable Velcro strap or an adjustable leather strap. The bottom of the sandal is slotted for ease in finding the proper mounting location. Sandals are available with either one or two straps and are sold in pairs. All mounting hardware is included. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.

B 3.3.1.- BioForm Sandals with one Strap w/Velcro



Model	Description	Size	Dimension	
			W	D
LIS	BioForm Sandals with one Strap w/Velcro, pair	Small	7	4
LIM	BioForm Sandals with one Strap w/Velcro, pair	Medium	9	4.5
LIL	BioForm Sandals with one Strap w/Velcro, pair	Large	11	5
LIX	BioForm Custom Sandals with one Strap w/Velcro	Custom	Template or Specs	

B 3.3.2.- BioForm Sandals with one Strap w/Buckle




Model	Description	Size	Dimension	
			W	D
LJS	BioForm Sandals with one Strap w/Buckle, pair	Small	7	4
LJM	BioForm Sandals with one Strap w/Buckle, pair	Medium	9	4.5
LJL	BioForm Sandals with one Strap w/Buckle, pair	Large	11	5
LJX	BioForm Custom Sandals with one Strap w/Buckle	Custom	Template or Specs	

“If you don’t see it... Contact Us”


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.3.3.- BioForm Sandals with one Leather Strap


Price	Model	Description	Size	Dimension	
				W	D
225	LMS	BioForm Sandals with one Leather Strap, pair	Small	7	4
225	LMM	BioForm Sandals with one Leather Strap, pair	Medium	9	4.5
225	LML	BioForm Sandals with one Leather Strap, pair	Large	11	5
Call	LMX	BioForm Custom Sandals with one Leather Strap	Custom	Template or Specs	

B 3.3.4.- BioForm Sandals with two Straps w/Velcro


Model	Description	Size	Dimension	
			W	D
LKS	BioForm Sandals with two Straps w/Velcro, pair	Small	7	4
LKM	BioForm Sandals with two Straps w/Velcro, pair	Medium	9	4.5
LKL	BioForm Sandals with two Straps w/Velcro, pair	Large	11	5
LKX	BioForm Custom Sandals with two Straps w/Velcro	Custom	Template or Specs	

B 3.3.5.- BioForm Sandals with two Straps w/Buckle


Model	Description	Size	Dimension	
			W	D
LLS	BioForm Sandals with two Straps w/Buckle, pair	Small	7	4
LLM	BioForm Sandals with two Straps w/Buckle, pair	Medium	9	4.5
LLL	BioForm Sandals with two Straps w/Buckle, pair	Large	11	5
LLX	BioForm Custom Sandals with two Straps w/Buckle	Custom	Template or Specs	


“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.3.6.- BioForm Sandals with two Leather Straps



Model	Description	Size	Dimension	
			W	D
LNS	BioForm Sandals with two Leather Straps, pair	Small	7	4
LNM	BioForm Sandals with two Leather Straps, pair	Medium	9	4.5
LNL	BioForm Sandals with two Leather Straps, pair	Large	11	5
LNX	BioForm Custom Sandals with two Leather Straps	Custom	Template or Specs	


B 3.4.- Solid Heel Cups

B 3.4.1.- BioForm Solid Heel Cups with Straps w/ Velcro

This choice for foot positioning is optimal when the primary goal is to secure the foot in a specific location on the footplate while allowing for some front of foot movement. Heel Cups are very beneficial for consumers with significant lower extremity tone.

BioForm Solid Heel Cups function in a similar fashion as the BioForm Sandals do. The base of the Heelcup is made from ½" HDPE, pre-drilled with countersunk holes for seamless mounting. The rear wall is made from ¼" kydex and is attached to the Heel Cup base with screws. The strap is affixed to the base of the Heelcup and is available with a Fastex Buckle closure or a Velcro closure.

Available in small, medium and large and come complete with all mounting hardware. Custom sizes available by selecting the appropriate model #. Sold in pairs.



BioForm Solid Heel Cups with Straps w/ Velcro


Model	Description	Size	Dimension		
			W	D	H
LAS	BioForm Solid Heel Cups with Straps w/ Velcro, pair	Small	4	7	2.5
LAM	BioForm Solid Heel Cups with Straps w/ Velcro, pair	Medium	4.5	9	2.5
LAL	BioForm Solid Heel Cups with Straps w/ Velcro, pair	Large	5	11	2.5
LAX	BioForm Custom Solid Heel Cups with Straps w/ Velcro	Custom	Template or Specs		

“If you don’t see it.... Contact Us”


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.4.2.- BioForm Solid Heel Cups with Straps w/ Fastex Buckle


Model	Description	Size	Dimension		
			W	D	H
LBS	BioForm Solid Heel Cups with Straps w/ Fastex Buckle, pair	Small	4	7	2.5
LBM	BioForm Solid Heel Cups with Straps w/ Fastex Buckle, pair	Medium	4.5	9	2.5
LBL	BioForm Solid Heel Cups with Straps w/ Fastex Buckle, pair	Large	5	11	2.5
LBX	BioForm Custom Solid Heel Cups with Straps w/ Fastex Buckle	Custom	Template or Specs		

B 3.5.- BioForm Foot Support Accessories
B 3.5.1.- BioForm Scuba Heel Loops


BioForm Scuba Heel Loops are an industry first. These “crush proof” heel loops are the answer to the chronic issue of heel loops getting crushed to the surface of a footplate. They are fabricated from a special webbing and come complete with all of the necessary mounting hardware. Available in small, medium and large and are sold in pairs. Custom sizes available by selecting the appropriate model #.


BioForm Scuba Heel Loops

Model	Description	Size	Dimension	
			L	H
LCS	BioForm Scuba Heel Loops, pair	Small	6	2
LCM	BioForm Scuba Heel Loops, pair	Medium	8	2
LCL	BioForm Scuba Heel Loops, pair	Large	10	2
LCX	BioForm Custom Scuba Heel Loops	Custom	Template or Specs	

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.5.2.- BioForm Ankle or Toe Straps

The BioForm Ankle or Toe Straps are an effective means of holding a foot or ankle in place upon a footplate. They are available in 2 styles: padded and non-padded. Also, the non-padded style comes with the choice of Velcro or Fastex Buckle closure. Each strap is 1” wide, has Tri-Glides for easy adjusting and comes with D-Rings that house the end fittings. Our end fittings provide a versatile tool for easy and convenient mounting. Available in small, medium and large and are sold in pairs. Custom sizes available by selecting the appropriate model #.

B 3.5.2.1.- BioForm Ankle or Toe Straps, Velcro Closure

	Model	Description	Size	Dimension	
				From	To
	MDS	BioForm Ankle or Toe Straps, 1”, D-Ring Velcro Closure, End fittings, pair	Small	7	10.5
	MDM	BioForm Ankle or Toe Straps, 1”, D-Ring Velcro Closure, End fittings, pair	Medium	8	13.5
	MDL	BioForm Ankle or Toe Straps, 1”, D-Ring Velcro Closure, End fittings, pair	Large	9	18.5
	MDX	BioForm Custom Ankle or Toe Straps, D-Ring Velcro Closure, End fittings	Custom	Template or Specs	

B 3.5.2.2.- BioForm Ankle or Toe Straps, Fastex Buckle Closure

	Model	Description	Size	Dimension	
				From	To
	MFS	BioForm Ankle or Toe Straps, 1”, Fastex Buckle Closure, End fittings, pair	Small	7	10.5
	MFM	BioForm Ankle or Toe Straps, 1”, Fastex Buckle Closure, End fittings, pair	Medium	8	13.5
	MFL	BioForm Ankle or Toe Straps, 1”, Fastex Buckle Closure, End fittings, pair	Large	9	18.5
	MFX	BioForm Custom Ankle or Toe Straps, Fastex Buckle Closure, End fittings	Custom	Template or Specs	

B 3.5.2.3.- BioForm Padded Ankle or Toe Straps, Fastex Buckle Closure

	Model	Description	Size	Dimension	
				From	To
	MPS-NB	BioForm Padded Ankle or Toe Straps, 1”, Fastex Buckle Closure, End fittings, pair	Small	7	10.5
	MPM-NB	BioForm Padded Ankle or Toe Straps, 1”, Fastex Buckle Closure, End fittings, pair	Medium	8	13.5
	MPL-NB	BioForm Padded Ankle or Toe Straps, 1”, Fastex Buckle Closure, End fittings, pair	Large	9	18.5
	MPX-NB	BioForm Custom Padded Ankle or Toe Straps, Fastex Buckle Closure, End fittings	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- MPS-NB becomes MPS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.5.2.4.- Pads for Ankle or Toe Straps

The BioForm Pads for Ankle or Toe Straps provide an opportunity to add padded protection for straps that were ordered without padding. The pads are sleeve type for easy installation.

BioForm Padded Sleeves for Ankle or Toe Straps

	<i>Model</i>	<i>Description</i>	<i>Dimension</i>	
			W	L
	LSP-NB	BioForm Padded Sleeves for Ankle or Toe Straps, pair	2	5.5
LSX-NB	BioForm Custom Padded Sleeves for Ankle or Toe Straps	Template or Specs		

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- LSP-NB becomes LSP-NE when color Eggplant is desired.

B 3.5.2.5.- BioForm Leather Straps

The BioForm Leather Straps are a versatile tool for very aggressive positioning needs of the foot, ankle or any other area with strong, unwanted movements. They are sold in pairs and are 1” in width. The straps are 2 pieces that are secured together with a Standard Duty Steel “belt type” buckle and have a hole at both ends for easy mounting.

BioForm Leather Straps

	Model	Description	Size	Dimension	
				From	To
	MMS	BioForm Leather Ankle or Toe Straps, 1” wide, pair	Small	12	16
	MMM	BioForm Leather Ankle or Toe Straps, 1” wide, pair	Medium	14	20
MML	BioForm Leather Ankle or Toe Straps, 1” wide, pair	Large	19	24	

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.- BioForm Footboxes



The BioForm line of Footboxes offer a wide array of styles and mounting options. Our 2 styles of Footboxes are: single leg boxes (designed to accommodate one leg per Footbox), and a 1 piece Footbox (designed to accommodate both legs in one Footbox). The Footboxes are available with hardware for: mounting to footplates, mounting to hanger tubes or mounting to a solid seat pan. We also offer our Footboxes without mounting hardware to offer you a more local/custom mounting solution. All Footbox options are available padded or non-padded. Padded units have the pad modules secured to the walls of the Footboxes with t-nuts and screws for added durability. Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#.

B 3.6.1.- BioForm, Double Leg (1 piece), Non-padded Footboxes

B 3.6.1.1.- BioForm, Double Leg (1 piece), Non-padded Footboxes, No Hardware





Model	Description	Dimension		
		W	D	H
EAS	BioForm Double Leg Footbox, No Hardware, for 14” Wheelchairs	11.5	8	8
EAM	BioForm Double Leg Footbox, No Hardware, for 16” Wheelchairs	13.5	9	9
EAL	BioForm Double Leg Footbox, No Hardware, for 18” Wheelchairs	15.5	10	10.5
EAO	BioForm Double Leg Footbox, No Hardware, for 20” Wheelchairs	17.5	11	12
EAX	BioForm Custom Double Leg Footbox, No Hardware	Template or Specs		

“If you don’t see it... Contact Us”



BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.1.2.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Nuts and Bolts to secure from footplates.


Model	Description	Dimension		
		W	D	H
EBS	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs	11.5	8	8
EBM	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs	13.5	9	9
EBL	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs	15.5	10	10.5
EBO	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs	17.5	11	12
EBX	BioForm Custom Double Leg Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

B 3.6.1.3.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Hardware to secure from hangers.




Model	Description	Dimension		
		W	D	H
ECS	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 14" Wheelchairs	11.5	8	8
ECM	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 16" Wheelchairs	13.5	9	9
ECL	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 18" Wheelchairs	15.5	10	10.5
ECO	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 20" Wheelchairs	17.5	11	12
ECX	BioForm Custom Double Leg Footbox, Hardware to secure from Hangers	Template or Specs		

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.1.4.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Hardware to secure from seat pan.


Model	Description	Dimension		
		W	D	H
EDS	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs	11.5	8	8
EDM	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs	13.5	9	9
EDL	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs	15.5	10	10.5
EDO	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs	17.5	11	12
EDX	BioForm Custom Double Leg Footbox, Hardware to secure from Seat Pan	Template or Specs		

B 3.6.2.- BioForm, Single Leg (2 pieces), Non-padded Footboxes
B 3.6.2.1.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, No Hardware



Model	Description	Dimension		
		W	D	H
EES	BioForm Single Leg Footbox, No Hardware, for 14" Wheelchairs (pair)	5.5	8	8
EEM	BioForm Single Leg Footbox, No Hardware, for 16" Wheelchairs (pair)	6.5	9	9
EEL	BioForm Single Leg Footbox, No Hardware, for 18" Wheelchairs (pair)	7.5	10	10.5
EEO	BioForm Single Leg Footbox, No Hardware, for 20" Wheelchairs (pair)	8.5	11	12
EEX	BioForm Custom Single Leg Footbox, No Hardware	Template or Specs		

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.2.2.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Nuts and Bolts to secure from footplates.


Model	Description	Dimension		
		W	D	H
EFS	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs (pair)	5.5	8	8
EFM	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs (pair)	6.5	9	9
EFL	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs (pair)	7.5	10	10.5
EFO	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs (pair)	8.5	11	12
EFX	BioForm Custom Single Leg Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

B 3.6.2.3.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Hardware to secure from hangers.




Model	Description	Dimension		
		W	D	H
EGS	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 14" Wheelchairs (pair)	5.5	8	8
EGM	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 16" Wheelchairs (pair)	6.5	9	9
EGL	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 18" Wheelchairs (pair)	7.5	10	10.5
EGO	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 20" Wheelchairs (pair)	8.5	11	12
EGX	BioForm Custom Single Leg Footbox, Hardware to secure from Hangers	Template or Specs		

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.2.4.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Hardware to secure from seat pan.


Model	Description	Dimension		
		W	D	H
EHS	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs (pair)	5.5	8	8
EHM	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs (pair)	6.5	9	9
EHL	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs (pair)	7.5	10	10.5
EHO	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs (pair)	8.5	11	12
EHX	BioForm Custom Single Leg Footbox, Hardware to secure from Seat Pan	Template or Specs		

B 3.6.3.- BioForm, Double Leg (1 piece), Padded Footboxes
B 3.6.3.1.- BioForm, Double Leg (1 piece), Padded Footboxes, No Hardware


Model	Description	Dimension		
		W	D	H
EIS-V393	BioForm Double Leg Padded Footbox, No Hardware, for 14" Wheelchairs	11.5	8	8
EIM-V393	BioForm Double Leg Padded Footbox, No Hardware, for 16" Wheelchairs	13.5	9	9
EIL-V393	BioForm Double Leg Padded Footbox, No Hardware, for 18" Wheelchairs	15.5	10	10.5
EIO-V393	BioForm Double Leg Padded Footbox, No Hardware, for 20" Wheelchairs	17.5	11	12
EIX-V393	BioForm Custom Double Leg Padded Footbox, No Hardware	Template or Specs		



If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- EIS-V393 becomes EIS-V354 when Vinyl Royal Blue is desired.

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.3.2.- BioForm, Double Leg (1 piece), Padded Footboxes, Nuts and Bolts to secure from footplates.


Model	Description	Dimension		
		W	D	H
EJS-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs	11.5	8	8
EJM-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs	13.5	9	9
EJL-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs	15.5	10	10.5
EJO-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs	17.5	11	12
EJX-V393	BioForm Custom Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- EJS-V393 becomes EJS-V354 when Vinyl Royal Blue is desired.

B 3.6.3.3.- BioForm, Double Leg (1 piece), Padded Footboxes, Hardware to secure from hangers.


Model	Description	Dimension		
		W	D	H
EKS-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 14" Wheelchairs	11.5	8	8
EKM-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 16" Wheelchairs	13.5	9	9
EKL-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 18" Wheelchairs	15.5	10	10.5
EKO-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 20" Wheelchairs	17.5	11	12
EKX-V393	BioForm Custom Double Leg Padded Footbox, Hardware to secure from Hangers	Template or Specs		


If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- EKS-V393 becomes EKS-V354 when Vinyl Royal Blue is desired.

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us



E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.3.4.- BioForm, Double Leg (1 piece), Padded Footboxes, Hardware to secure from seat pan.


Model	Description	Dimension		
		W	D	H
ELS-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs	11.5	8	8
ELM-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs	13.5	9	9
ELL-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs	15.5	10	10.5
ELO-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs	17.5	11	12
ELX-V393	BioForm Custom Double Leg Padded Footbox, Hardware to secure from Seat Pan	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- ELS-V393 becomes ELS-V354 when Vinyl Royal Blue is desired.

B 3.6.4.- BioForm, Single Leg (2 pieces), Padded Footboxes
B 3.6.4.1.- BioForm, Single Leg (2 pieces), Padded Footboxes, No Hardware


Model	Description	Dimension		
		W	D	H
EMS-V393	BioForm Single Leg Padded Footbox, No Hardware, for 14" Wheelchairs (pair)	5.5	8	8
EMM-V393	BioForm Single Leg Padded Footbox, No Hardware, for 16" Wheelchairs (pair)	6.5	9	9
EML-V393	BioForm Single Leg Padded Footbox, No Hardware, for 18" Wheelchairs (pair)	7.5	10	10.5
EMO-V393	BioForm Single Leg Padded Footbox, No Hardware, for 20" Wheelchairs (pair)	8.5	11	12
EMX-V393	BioForm Custom Single Leg Padded Footbox, No Hardware	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- EMS-V393 becomes EMS-V354 when Vinyl Royal Blue is desired.

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.4.2.- BioForm, Single Leg (2 pieces), Padded Footboxes, Nuts and Bolts to secure from footplates.


Model	Description	Dimension		
		W	D	H
ENS-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs (pair)	5.5	8	8
ENM-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs (pair)	6.5	9	9
ENL-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs (pair)	7.5	10	10.5
ENO-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs (pair)	8.5	11	12
ENX-V393	BioForm Custom Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- ENS-V393 becomes ENS-V354 when Vinyl Royal Blue is desired.

B 3.6.4.3.- BioForm, Single Leg (2 pieces), Padded Footboxes, Hardware to secure from hangers.


Model	Description	Dimension		
		W	D	H
EOS-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 14" Wheelchairs (pair)	5.5	8	8
EOM-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 16" Wheelchairs (pair)	6.5	9	9
EOL-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 18" Wheelchairs (pair)	7.5	10	10.5
EOO-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 20" Wheelchairs (pair)	8.5	11	12
EOX-V393	BioForm Custom Single Leg Padded Footbox, Hardware to secure from Hangers	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- EOS-V393 becomes EOS-V354 when Vinyl Royal Blue is desired.



“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.6.4.4.- BioForm, Single Leg (2 pieces), Padded Footboxes, Hardware to secure from seat pan.



Model	Description	Dimension		
		W	D	H
EPS-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs (pair)	5.5	8	8
EPM-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs (pair)	6.5	9	9
EPL-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs (pair)	7.5	10	10.5
EPO-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs (pair)	8.5	11	12
EPX-V393	BioForm Custom Single Leg Padded Footbox, Hardware to secure from Seat Pan	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- EPS-V393 becomes EPS-V354 when Vinyl Royal Blue is desired.

"If you don't see it.... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.7.- BioForm Calf Support

B 3.7.1.- BioForm Soft Calf Panels

BioForm Soft Calf Panels are a useful option when assistance is needed in keeping the consumers' legs from sliding rearwards off of the footplates.

The line of BioForm Soft Calf Panels are available in 2 styles: Non-Stretch, Vinyl (Naugahyde) for a more firm support and Stretch, double layer of Neoprene for a panel that is softer and has more flexibility. Each comes with straps for easy securing to the footrest/legrest assembly and close with Velcro or a Fastex Buckle.

Available in small, medium, large and x-large. Custom sizes available by selecting the appropriate model #.

Each style is available in 4 sizes.

B 3.7.1.1.- BioForm Stretch Soft Black Calf Panels with Fastex Buckle Closure



Model	Description	Size	Dimension	
			W	H
NCS-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	Small	16	7
NCM-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	Medium	18	8
NCL-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	Large	20	8
NCO-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	X-Large	22	9
NCX-NB	BioForm Custom Stretch Soft Black Calf Panel, Layers of Neoprene, Fastex Buckles	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- NCS-NB becomes NCS-NE when color Eggplant is desired.

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.7.1.2.- BioForm Stretch Soft Black Calf Panels with Velcro Closure (4+2 straps)


Model	Description	Size	Dimension	
			W	H
NDS-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	Small	16	7
NDM-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	Medium	18	8
NDL-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	Large	20	8
NDO-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	X-Large	22	9
NDX-NB	BioForm Custom Stretch Soft Black Calf Panel, Layers of Neoprene, Straps w/ Velcro Closure	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- NDS-NB becomes NDS-NE when color Eggplant is desired.

B 3.7.1.3.- BioForm Non-Stretch Soft Black Calf Panels with Fastex Buckle Closure


Model	Description	Size	Dimension	
			W	H
VCS-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	Small	16	7
VCN-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	Medium	18	8
VCL-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	Large	20	8
VCO-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	X-Large	22	9
VCX-V393	BioForm Custom Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), Straps w/ Fastex Buckles	Custom	Template or Specs	

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- VCS- V393 becomes VCS-V322 when color Sun Yellow is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.7.1.4.- BioForm Non-Stretch Soft Black Calf Panels with Velcro Closure (4+2 straps)


Model	Description	Size	Dimension	
			W	H
VDS-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	Small	16	7
VDM-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	Medium	18	8
VDL-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	Large	20	8
VDO-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	X-Large	22	9
VDX-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), Straps w/ Velcro Closure	Custom	Template or Specs	

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- VDS- V393 becomes VDS-V322 when color Sun Yellow is desired.

B 3.7.2.- BioForm Solid Calf Panels

The BioForm Calf Panels are constructed from a ½” HDPE base, ½” Sunmate Foam, cover and mounting hardware. They are available with fixed mounting hardware as well as removable hardware.

B 3.7.2.1.- BioForm Solid Calf Panel, Fixed Mounting Hardware


Model	Description	Size	Dimension	
			W	H
CPS-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 13” & 14” Wheelchairs	Small	16	7
CPM-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 15” & 16” Wheelchairs	Medium	18	8
CPL-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 17” & 18” Wheelchairs	Large	20	8
CPO-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 19” & 20” Wheelchairs	X-Large	22	9
CPX-NB	BioForm Custom Solid Calf panel, fixed hardware	Custom	Template or Specs	


If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- CPS-NB becomes CPS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.7.2.2.- BioForm Solid Calf Panel, Removable Hardware


Model	Description	Size	Dimension	
			W	H
CRS-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, removable hardware, for 13” & 14” Wheelchairs	Small	16	7
CRM-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, removable hardware, for 15” & 16” Wheelchairs	Medium	18	8
CRL-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, removable hardware, for 17” & 18” Wheelchairs	Large	20	8
CRO-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, removable hardware, for 19” & 20” Wheelchairs	X-Large	22	9
CRX-NB	BioForm Custom Solid Calf panel, removable hardware	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.- CRS-NB becomes CRS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 3.8.- BioForm Anterior Knee Blocks

BioDynamics' BioForm Anterior Knee Blocks are designed to help properly position the lower extremities while assisting in stabilizing the pelvis. They are available in 2 styles of mounting options; center mount (removable) and center and lateral mount (removable with lateral straps for added security).

Anterior Knee Blocks are not recommended for clients with hip or knee problems and should not be utilized to hold the pelvis to the rear of the wheelchair.

Available in small, medium and large. Custom sizes available by selecting the appropriate model #.



Anterior Knee Blocks

Model	Description	Size	Dimension			
			W	H	WK	DK
CKS-NB	BioForm Anterior Knee Block, Center Mount, 14" wide.	Small	14	4	3	1
CKM-NB	BioForm Anterior Knee Block, Center Mount, 16" wide.	Medium	16	4.5	4	1
CKL-NB	BioForm Anterior Knee Block, Center Mount, 18" wide.	Large	18	5	5	1.5
CLS-NB	BioForm Anterior Knee Block, Center & Lateral Mount, 14" wide.	Small	14	4	3	1
CLM-NB	BioForm Anterior Knee Block, Center & Lateral Mount, 16" wide.	Medium	16	4.5	4	1
CLL-NB	BioForm Anterior Knee Block, Center & Lateral Mount, 18" wide.	Large	18	5	5	1.5
CKX-NB	BioForm Custom Anterior Knee Block to Spec.,	Custom	Template or Specs			

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- CKS-NB becomes CKS-NE when color Eggplant is desired.

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 4.- Headrest Systems



B 4.1.- Headrest Pads

BioDynamics manufactures a line of headrests which are a true complement to the line of seating and positioning systems. The line of headrests comes standard with a 1” layer of Sunmate foam and a covering material. Each material has properties which will make the selection unique for each clinical need (I.E.- constant friction, heat build up, etc.).

Additional covers are also available. Custom modifications to the headrests are available such as occipital contouring, ear relief cutouts, etc...

Flat Headrest Pads are an effective means of protecting the consumer’s head/neck during bus/ambulette transportation. Typically ordered with the removable extension bracket, these pads can be easily removed while in school or at work and quickly attached for transportation. A headrest hardware option needs to be ordered with this pad in order for it to be mounted to a back support.

BioCurve Headrest Pads are a contoured pad suitable for the client whose headrest needs will not be satisfied by the Flat Headrest Pad. This headrest’s function is enhanced with an optional occipital contour which will help maintain a neutral neck while supporting the head. A headrest hardware option needs to be ordered with this pad in order for it to be mounted to a back support.

BioForm Cervical Control Headrest Pads are similar in support to the BioCurve pads with a cervical control component. The ends of the pads have a “tail” which reaches down from the head area towards the neck offering cervical control. A headrest hardware option needs to be ordered with this pad in order for it to be mounted to a back support.



BioTri-Pad Headrest Pads have the ability to offer differing levels of control on either side of the head. This headrest selection is offered with angle adjustable head wings or fixed head wings. With the angle adjustable head wings, the care taker can move one wing in towards midline more than the other. The fixed angle head wings offer a cost effective option to procure a Tri-Pad headrest when significant funding restrictions exist. This would be the best option for a head that tends to lean/roll to one side. A headrest hardware option needs to be ordered with this pad in order for it to be mounted to a back support.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 4.1.1.- Flat Headrest Pad



Model	Description	Size	Dimension	
			W	H
HFS-NB	Flat Headrest Pad	Small	8	5
HFM- NB	Flat Headrest Pad	Medium	10	6
HFL- NB	Flat Headrest Pad	Large	12	7
HFO- NB	Flat Headrest Pad	X-Large	14	8
HFX- NB	Flat Headrest Pad	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.- HFS-NB becomes HFS-NE when color Eggplant is desired.

B 4.1.2.- BioCurve Headrest Pad


Model	Description	Size	Dimension	
			W	H
HWM- NB	BioCurve Headrest Pad	Medium	10	5.5
HWL- NB	BioCurve Headrest Pad	Large	13.5	6
HWX- NB	BioCurve Custom Headrest Pad	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.- HWM-NB becomes HWM-NE when color Eggplant is desired.

B 4.1.3.- BioForm Cervical Control Headrest Pad


Model	Description	Size	Dimension		
			A	B	C
HRS- NB	BioForm Cervical Control Headrest Pad	Small	7	5.5	12
HRM- NB	BioForm Cervical Control Headrest Pad	Medium	7	5.5	13.5
HRL- NB	BioForm Cervical Control Headrest Pad	Large	7.5	5.5	15.5
HRX- NB	BioForm Custom Cervical Control Headrest Pad	Custom	Template or Specs		


If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.- HRS-NB becomes HRS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 4.1.4.- BioForm Tri-Pad Headrest with Angle Adjustable Wings


Model	Description	Size	Dimension (ea. pc.)	
			W	H
HTS- NB	BioForm Tri-Pad Headrest with Angle Adjustable Wings	Small	5	5
HTM- NB	BioForm Tri-Pad Headrest with Angle Adjustable Wings	Medium	6	6
HTL- NB	BioForm Tri-Pad Headrest with Angle Adjustable Wings	Large	7	7
HTO- NB	BioForm Tri-Pad Headrest with Angle Adjustable Wings	X-Large	8	8
HTX- NB	BioForm Custom Tri-Pad Headrest with Angle Adjustable Wings	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- HTS-NB becomes HTS-NE when color Eggplant is desired.

B 4.1.5.- BioForm Tri-Pad Headrest with Fixed Angle Wings


Model	Description	Size	Dimension (ea. pc.)	
			W	W
HGA-NB	BioForm Tri-Pad Headrest with Fixed Angle Wings	Small	5	5
HGB-NB	BioForm Tri-Pad Headrest with Fixed Angle Wings	Medium	6	6
HGC-NB	BioForm Tri-Pad Headrest with Fixed Angle Wings	Large	7	7
HGX-NB	BioForm Custom Tri-Pad Headrest with Fixed Angle Wings	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- HTS-NB becomes HTS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 4.1.6.- BioCurve Headrest Package



HWMIC



HWLIC

Model	Description	Size	Dimension	
			W	W
HWMIC- NB	BioCurve Headrest Pad with removable Height, Depth and Angle Adjustable Hardware	Medium	10	5.5
HWLIC- NB	BioCurve Headrest Pad with removable Height, Depth and Angle Adjustable Hardware	Large	13.5	6

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.-HWMIC-NB becomes HWMIC-NE when color Eggplant is desired.

B 4.2.- Headrest Support Modifications

Model	Description
HMC	Occipital Contour
HME	Ear Relief Cut-Out
HMR-NB	Removable Cover (1 st)
HMS-NB	Extra Removable Cover
HMX	Custom Modification

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.
i.e.-HMR-NB becomes HMR-NE when color Eggplant is desired.

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 4.3.- Headrest Support Mounting Hardware

BioDynamics offers a variety of mounting options for headrest pads.

- Removable Extension Hardware offers vertical height adjustment of the pad selected. This bracket will allow the headrest pad to be quickly and easily removed by loosening a knob.
- Removable Height and Depth Adjustable Hardware is made from square, stainless steel stock. The choice of using stainless steel is to prevent rusting as well as offering maximum strength/durability. The square shape prevents rotating of the horizontal and vertical bars. The height and depth of the vertical and horizontal rods, respectively, can be adjusted. The angle of the pad relative to the bracket is fixed. This bracket will allow the headrest pad to be quickly and easily removed by loosening a knob.
- Removable Height, Depth and Angle Adjustable Headrest Hardware is made from round, stainless steel stock. The stainless steel material will prevent rusting and offer maximum durability. The height and depth of the vertical and horizontal rods can be adjusted as well as the angle of the pad relative to the bracket. This bracket also allows the headrest pad to be quickly and easily removed by loosening a knob.
- Flip-Down Height, Depth and Angle Adjustable Headrest Hardware is made from round, stainless steel stock. The stainless steel material will prevent rusting and offer maximum durability. The height and depth of the vertical and horizontal rods can be adjusted as well as the angle of the pad relative to the bracket. This bracket flips down with the pull of a spring loaded pin. Selecting this option gives the caretaker the benefit of moving the headrest pad out of the way for transfers, etc., while keeping it attached to the chair. This is beneficial for loss prevention.

Headrest Support Mounting Hardware



#HIB



HID



#HIC



#HIE

Model	Description
HIB	Removable Extension Headrest Hardware
HID	Removable Height and Depth Adjustable Hardware
HIC	Removable Height, Depth and Angle Adjustable Hardware
HIE	Flip-Down, Height, Depth and Angle Adjustable Hardware
HIX	Custom Hardware

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 5.- Upper Extremity Supports

Upper Extremity Supports offered by Bio Dynamics are designed to enhance the activities of daily living (adl's) as well as aid in upper trunk support. The many modifications available will help create the final product that will suit all of the client's needs.

Upper Extremity Support manufactured by BioDynamics are offered in various forms; arm troughs, upper extremity support surfaces, arm/wrist straps, etc.,

B 5.1.- Upper Extremity Support Surface



Upper Extremity Supports Surfaces offered by Bio Dynamics are designed to enhance the activities of daily living (adl's) as well as aid in upper trunk support. The many modifications available will help create the final product that will suit all of the client's needs. The various choices of mounting hardware allow our support surfaces to fit on most any mobility base.

The Upper Extremity Supports Surfaces are available in 3 surface materials; 3/8" clear Lexan, 1/4" Black Phenolic, 3/8" White Sanalite and 3/8" Black Sanalite.

BioDynamics offers 2 lines of Upper Extremity Supports Surfaces: BioDezine and BioForm. Four packages exist within both of the: BioDezine and BioForm lines of upper extremity support surfaces.

The list of modifications available to our line of Upper Extremity Supports Surfaces is detailed below.

The dimensions of the Upper Extremity Supports Surfaces need to be specified.

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 5.1.1.- BioDezine Upper Extremity Support Surface, Package

The BioDezine Upper Extremity Supports Surfaces are available in 4 surface materials; 3/8” clear Lexan, 1/4” Black Phenolic, 3/8” White Sanalite and 3/8” Black Sanalite. Our standard mounting option of slides with an attaching belt comes with an adjustable stop which controls how close the inner rim of the support surface will come to the client. The BioDezine Upper Extremity Support Package is also offered standard with up to 3 rims on the perimeter of the surface.





For Custom Trays



For Standard Trays

BioDezine Upper Extremity Support Surface, Package

Model	Description	Size	Dimension				
			W	D	WC	DC	
 CAA to CAD	CAA	BioDezine Upper Extremity Support Tray, 3/8” Clear Lexan, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Small	20	18	10	7
	CAB	BioDezine Upper Extremity Support Tray, 3/8” Clear Lexan, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Medium	22	20	12	9
	CAC	BioDezine Upper Extremity Support Tray, 3/8” Clear Lexan, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Large	24	24	16	11
	CAD	Custom BioDezine Upper Extremity Support Tray, 3/8” Clear Lexan, up to 24”x24”, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Custom	Template or Specs			
 CAE to CAH	CAE	BioDezine Upper Extremity Support Tray, 1/4” Black Phenolic, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Small	20	18	10	7
	CAF	BioDezine Upper Extremity Support Tray, 1/4” Black Phenolic, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Medium	22	20	12	9
	CAG	BioDezine Upper Extremity Support Tray, 1/4” Black Phenolic, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Large	24	24	16	11
	CAH	Custom BioDezine Upper Extremity Support Tray, 1/4” Black Phenolic, up to 24”x24”, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Custom	Template or Specs			
 CAI to CAL	CAI	BioDezine Upper Extremity Support Tray, 3/8” White Sanalite, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Small	20	18	10	7
	CAJ	BioDezine Upper Extremity Support Tray, 3/8” White Sanalite, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Medium	22	20	12	9
	CAK	BioDezine Upper Extremity Support Tray, 3/8” White Sanalite, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Large	24	24	16	11
	CAL	Custom BioDezine Upper Extremity Support Tray, 3/8” White Sanalite, up to 24”x24”, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Custom	Template or Specs			

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


 <p>CAM to CAP</p>	CAM	BioDezine Upper Extremity Support Tray, 3/8" Black Sanalite, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Small	20	18	10	7
	CAN	BioDezine Upper Extremity Support Tray, 3/8" Black Sanalite, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Medium	22	20	12	9
	CAO	BioDezine Upper Extremity Support Tray, 3/8" Black Sanalite, Slides w/ adjustable tray stops, Belt & up to 3 Rims	Large	24	24	16	11
	CAP	Custom BioDezine Upper Extremity Support Tray, 3/8" Black Sanalite, up to 24"x24", Slides w/ adjustable tray stops, Belt & up to 3 Rims	Custom	Template or Specs			

B 5.1.2.- BioForm Upper Extremity Support Surface, Package

The BioForm Upper Extremity Supports Surfaces are available in 4 surface materials; 3/8" clear Lexan, 1/4" Black Phenolic, 3/8" White Sanalite and 3/8" Black Sanalite. Our standard mounting option of slides with an attaching belt insures the Upper Extremity Support Surface will stay in place.



For Custom Trays



For Standard Trays

BioForm Upper Extremity Support Surface, Package

Model	Description	Size	Dimension				
			W	D	WC	DC	
 <p>CBA to CBD</p>	CBA	BioForm Upper Extremity Support Tray, 3/8" Clear Lexan, Slides & Belt	Small	20	18	10	7
	CBB	BioForm Upper Extremity Support Tray, 3/8" Clear Lexan, Slides & Belt	Medium	22	20	12	9
	CBC	BioForm Upper Extremity Support Tray, 3/8" Clear Lexan, Slides & Belt	Large	24	24	16	11
	CBD	Custom BioForm Upper Extremity Support Tray, 3/8" Clear Lexan, up to 24"x24", Slides & Belt	Custom	Template or Specs			
 <p>CBE to CBH</p>	CBE	BioForm Upper Extremity Support Tray, 1/4" Black Phenolic, Slides & Belt	Small	20	18	10	7
	CBF	BioForm Upper Extremity Support Tray, 1/4" Black Phenolic, Slides & Belt	Medium	22	20	12	9
	CBG	BioForm Upper Extremity Support Tray, 1/4" Black Phenolic, Slides & Belt	Large	24	24	16	11
	CBH	Custom BioForm Upper Extremity Support Tray, 1/4" Black Phenolic, up to 24"x24", Slides & Belt	Custom	Template or Specs			

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

 CBI to CBL	CBI	BioForm Upper Extremity Support Tray, 3/8" White Sanalite, Slides & Belt	Small	20	18	10	7
	CBJ	BioForm Upper Extremity Support Tray, 3/8" White Sanalite, Slides & Belt	Medium	22	20	12	9
	CBK	BioForm Upper Extremity Support Tray, 3/8" White Sanalite, Slides & Belt	Large	24	24	16	11
	CBL	Custom BioForm Upper Extremity Support Tray, 3/8" White Sanalite, up to 24"x24", Slides & Belt	Custom	Template or Specs			
 CBM to CBP	CBM	BioForm Upper Extremity Support Tray, 3/8" Black Sanalite, Slides & Belt	Small	20	18	10	7
	CBN	BioForm Upper Extremity Support Tray, 3/8" Black Sanalite, Slides & Belt	Medium	22	20	12	9
	CBO	BioForm Upper Extremity Support Tray, 3/8" Black Sanalite, Slides & Belt	Large	24	24	16	11
	CBP	Custom BioForm Upper Extremity Support Tray, 3/8" Black Sanalite, up to 24"x24", Slides & Belt	Custom	Template or Specs			


B 5.1.3.- Upper Extremity Support Surface Materials

The BioDynamics Upper Extremity Supports Surface materials have been selected to maximize the life span and integrity of the support surface.

Clear Lexan offers the consumer the ability to enjoy the benefits of being able to view their lower extremities while offering a very durable surface.

The Black Phenolic surface is a pulp based product that offers a smooth, lightweight, non transparent surface.

The Black and White Sanalite is an FDA approved surface for the food service industry making it an ideal selection for "messy" applications and/or situations which require an easy to clean surface.



Upper Extremity Support Surface Materials

Model	Description	Dimension
CUA	3/8" Clear Lexan, up to 24" x 24"	Template or Specs
CUB	1/4" Black Phenolic, up to 24" x 24"	Template or Specs
CUC	3/8" White Sanalite, up to 24" x 24"	Template or Specs
CUD	3/8" Black Sanalite, up to 24" x 24"	Template or Specs

"If you don't see it... Contact Us"


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 5.1.4.- Upper Extremity Supports Surface Modifications

The extensive offering of modifications listed below can be selected in conjunction with any of the Upper Extremity Support Surface Packages or Upper Extremity Support Surface Materials.



Upper Extremity Supports Surface Modifications

Model	Description	Dimension
CUF	Custom Perimeter Cut-outs (I.E.-Joystick, Power seat function switch, self propelling notches, etc.)	Template or Specs
CUG	Custom Surface Cut-outs (Joystick, Power seat function switch, etc.)	Template or Specs
CUH	Oversize width (from 24" to 30")	Template or Specs
CUI	Oversize depth (from 24" to 30")	Template or Specs
CUJ	Set of Rims (up to 3)	Template or Specs
CUK-V393	Pad elbow portion of board *	Template or Specs
CUV-V393	Pad Inner rim of board *	Template or Specs
CUM-V393	Pad Top of board *	Template or Specs
CUN-V393	Pad Bottom of board *	Template or Specs
CUO-V393	Padded Protraction Blocks w/Hardware (pair) *	Template or Specs
CUP	Clear Plastic Overlay with Velcro attachment	Template or Specs
CUU	Clear Plastic Underlay with Velcro attachment	Template or Specs
CUW	Upgrade Belt to Auto Buckle Style	1.5"
CUX	Custom Upper Extremity Support Modification	Template or Specs

Note:

* All padded modifications are covered with **black vinyl** unless otherwise requested

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- CUK-V393 becomes CUK-V354 when Vinyl Royal Blue is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 5.1.5.- Upper Extremity Support Surface Mounting Hardware

Our standard mounting option of slides with an attaching belt comes with an adjustable stop which controls how close the inner rim of the support surface will come to the client.

The other option is a single slide on one side and a limiting wall with a stop and two sets of velcro straps on the other side. This option is useful when used with a power wheelchair. It allows for the support surface to be placed on to the arms instead of sliding which would be hindered by the joystick.



Upper Extremity Support Surface Mounting Hardware

Model	Description	Dimension
CHA	Two Slides and a Belt with Fastex Buckle	Template or Specs
CHB	Two Slides with built-in adjustable stops and a Belt with Fastex Buckle	Template or Specs
CHR	Single Slide and a limiting wall on Rt side	Template or Specs
CHL	Single Slide and a limiting wall on Lf side	Template or Specs

“If you don’t see it.... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


B 6.- Miscellaneous Accessories

Miscellaneous Accessories consist of a variety of aftermarket handles to fit on most any power wheelchair's joystick. Each selection enables users of different abilities access to power mobility.

B 6.1.- BioForm Joystick Handles

B 6.1.1.- BioForm Goal Post Joystick Handles

BioForm Goal Post joystick handles ('Y' shape), provides a superb option for greater control of a power wheelchair for the individual who lacks dexterity with their hands. They are shaped from a piece of 1/16" x 1" aluminum with a coupling piece on the underside for easy attaching to the joystick. This coupling device has a set screw for securing it to the joystick.




BioForm Goal Post Joystick Handle

Model	Description	Size	Dimension		
			W	GH	H
JGS	Goal Post Joystick Handle w/set screw	Small	3	1	2.5
JGM	Goal Post Joystick Handle w/set screw	Medium	4	1	3
JGL	Goal Post Joystick Handle w/set screw	Large	5	1	3.5
JGX	Custom Goal Post Joystick Handle	Custom	Template or Specs		

B 6.1.2.- BioForm 'T' Joystick Handles

BioForm 'T' joystick handles provide an option for excellent control of a power wheelchair for the individual who lacks dexterity with their hands. They are shaped from a 1/2" aluminum rod with a coupling piece on the underside for easy attaching to the joystick. This coupling device has a set screw for securing it to the joystick.



BioForm 'T' Joystick Handle

Model	Description	Size	Dimension		
			W	D	H
JTS	'T' Joystick Handle w/set screw, 3" wide	Small	3	0.5	2.5
JTM	'T' Joystick Handle w/set screw, 4" wide	Medium	4	0.5	3
JTL	'T' Joystick Handle w/set screw, 5" wide	Large	5	0.5	3.5
JTX	Custom 'T' Joystick Handle	Custom	Template or Specs		

"If you don't see it.... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.