

B 2.- Positioning Belts, Straps and Pads

Positioning Belts, Straps and Pads found in this section enhance the positioning of the upper trunk and lower extremities. Many of the Positioning Belts and Straps are available with pads. Those without can be retrofitted with our padded sleeve options. All choices come with the appropriate mounting hardware.

B 2.1- Anterior Trunk Supports



The BioForm Line of Anterior Trunk Supports are available in Neoprene (Rubatex), Reverse Dartex, Vinyl (Naugahyde), and Ballistic Nylon. Neoprene is effective when control is needed but a stretch component would prove beneficial. Reverse Dartex will offer a more rigid support than Neoprene while still having some elasticity, as well as providing an easy to clean surface. Vinyl (Naugahyde) will offer the most rigid support as well as a surface which handles fluids and cleanings most easily. Ballistic Nylon is the preferred material when the consumer will try to orally or otherwise try to destroy the harness.

All styles below are available in standard sizes from XS through XL. A customization to the size and/or shape is available by selecting the appropriate model #.

“If you don’t see it... Contact Us”


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

- **BioForm Standard Anterior Trunk Supports** are an effective way to help maintain a properly positioned trunk when significant assistance is needed. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm Anterior Trunk Supports with Snap Closures** assist in easier/quicker transfers from the wheelchair. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm “H” type Anterior Trunk Supports** are effective when less than maximum assistance is needed to maintain upper trunk position and quick and easy transfer is a primary goal. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm Low Profile Anterior Trunk Supports** are effective when significant assistance is needed to properly maintain upper trunk position and a g-tube needs to be accommodated for. This product may also prove helpful when used with the adolescent to fully matured female shape. Comes complete with 4 securing straps, end fittings and all mounting hardware.
- **BioForm 6 Point Trunk Support** is effective for the client who needs **maximum** trunk support: Anterior Trunk, pelvis and shoulders. Comes complete with 6 securing straps, end fittings and all mounting hardware.

B 2.1.1.- BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support



Model	Description	Size	Dimension					
			A	B	C	D	E	F
NRT-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
NRD-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
NRS-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	Small	17	9	10	5	2	5
NRM-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
NRL-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
NRO-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
NRX-NB	BioForm Stretch (Black Neoprene) Custom Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NRD-NB becomes NRD-NE when color Eggplant is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.2.- BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
NLD-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	X-Small	14	8	8.5	4	2	1
NLS-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	Small	17	9	10	5	2	1
NLM-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	Medium	20	10	11.5	6	2	1.25
NLL-NB	BioForm Stretch (Black Neoprene) Low-Profile Anterior Trunk Support	Large	22	11	12	7	2.5	1.5
NLX-NB	BioForm Stretch (Black Neoprene) Custom Low-Profile Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NLD-NB becomes NLD-NE when color Eggplant is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.3.- BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure


Model	Description	Size	Dimension					
			A	B	C	D	E	F
NST-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	XX-Small	10	5	5	4	1.75	1.25
NSD-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	X-Small	14	8	8.5	4	2	4
NSS-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	Small	17	9	10	5	2	5
NSM-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	Medium	20	10	11.5	6	2	6
NSL-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	Large	22	11	12	7	2.5	7
NSO-NB	BioForm Stretch (Black Neoprene) Standard Anterior Trunk Support, Snap Closure	X-Large	24	12	14	8	2.5	8
NSX-NB	BioForm Stretch (Black Neoprene) Custom Anterior Trunk Support, Snap Closure	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NSD-NB becomes NSD-NE when color Eggplant is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.4.- BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support


Model	Description	Size	Dimension				
			A	B	C	D	E
NHS-NB	BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support	Small	10	14	6	8	2
NHM-NB	BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support	Medium	13	18	8	10	2.25
NHL-NB	BioForm Stretch (Black Neoprene) “H”-Style Anterior Trunk Support	Large	16	21	10	12	2.5
NHX-NB	BioForm Stretch (Black Neoprene) Custom “H”-Style Anterior Trunk Support	Custom	Template or Specs				

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NHD-NB becomes NHD-NE when color Eggplant is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.5.- BioForm Stretch (Black Neoprene) 6 Point Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
NMT-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	XX-Small	10	5	5	4	1.75	1.25
NMD-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	X-Small	14	8	8.5	4	2	4
NMS-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	Small	17	9	10	5	2	5
NMM-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	Medium	20	10	11.5	6	2	6
NML-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	Large	22	11	12	7	2.5	7
NMO-NB	BioForm Stretch (Black Neoprene) 6 Point Trunk Support	X-Large	24	12	14	8	2.5	8
NMX-NB	BioForm Stretch (Black Neoprene) Custom 6 Point Trunk Support	Custom	Template or Specs					

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.
i.e.- NMD-NB becomes NMD-NE when color Eggplant is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.6.- BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
DRT	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
DRD	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
DRS	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	Small	17	9	10	5	2	5
DRM	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
DRL	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
DRO	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
DRX	BioForm Semi-Stretch (Black Reverse Dartex) Custom Anterior Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.7.- BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support



Model	Description	Size	Dimension					
			A	B	C	D	E	F
DLD	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	X-Small	14	8	8.5	4	2	1
DLS	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	Small	17	9	10	5	2	1
DLM	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	Medium	20	10	11.5	6	2	1.25
DLL	BioForm Semi-Stretch (Black Reverse Dartex) Low-Profile Anterior Trunk Support	Large	22	11	12	7	2.5	1.5
DLX	BioForm Semi-Stretch (Black Reverse Dartex) Custom Low-Profile Anterior Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

B 2.1.8.- BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure


Model	Description	Size	Dimension					
			A	B	C	D	E	F
DST	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	XX-Small	10	5	5	4	1.75	1.25
DSD	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	X-Small	14	8	8.5	4	2	4
DSS	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	Small	17	9	10	5	2	5
DSM	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	Medium	20	10	11.5	6	2	6
DSL	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	Large	22	11	12	7	2.5	7
DSO	BioForm Semi-Stretch (Black Reverse Dartex) Standard Anterior Trunk Support, Snap Closure	X-Large	24	12	14	8	2.5	8
DSX	BioForm Semi-Stretch (Black Reverse Dartex) Custom Anterior Trunk Support, Snap Closure	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.9.- BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support


Model	Description	Size	Dimension				
			A	B	C	D	E
DHS	BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support	Small	10	14	6	8	2
DHM	BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support	Medium	13	18	8	10	2.25
DHL	BioForm Semi-Stretch (Black Reverse Dartex) “H”-Style Anterior Trunk Support	Large	16	21	10	12	2.5
DHX	BioForm Semi-Stretch (Black Reverse Dartex) Custom “H”-Style Anterior Trunk Support	Custom	Template or Specs				


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us


E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.10.- BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
DMT	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	XX-Small	10	5	5	4	1.75	1.25
DMD	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	X-Small	14	8	8.5	4	2	4
DMS	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	Small	17	9	10	5	2	5
DMM	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	Medium	20	10	11.5	6	2	6
DML	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	Large	22	11	12	7	2.5	7
DMO	BioForm Semi-Stretch (Black Reverse Dartex) 6 Point Trunk Support	X-Large	24	12	14	8	2.5	8
DMX	BioForm Semi-Stretch (Black Reverse Dartex) Custom 6 Point Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.11.- BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VRT-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
VRD-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
VRS-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	Small	17	9	10	5	2	5
VRM-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
VRL-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
VRO-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
VRX-V393	BioForm Non-Stretch (Black Vinyl) Custom Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VRD-V393 becomes VRD-V322 when color Sun Yellow is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.12.- BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VLD-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	X-Small	14	8	8.5	4	2	1
VLS-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	Small	17	9	10	5	2	1
VLM-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	Medium	20	10	11.5	6	2	1.25
VLL-V393	BioForm Non-Stretch (Black Vinyl) Low-Profile Anterior Trunk Support	Large	22	11	12	7	2.5	1.5
VLX-V393	BioForm Non-Stretch (Black Vinyl) Custom Low-Profile Anterior Trunk Support	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VLD-V393 becomes VLD-V322 when color Sun Yellow is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.13.- BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VST-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	XX-Small	10	5	5	4	1.75	1.25
VSD-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	X-Small	14	8	8.5	4	2	4
VSS-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	Small	17	9	10	5	2	5
VSM-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	Medium	20	10	11.5	6	2	6
VSL-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	Large	22	11	12	7	2.5	7
VSO-V393	BioForm Non-Stretch (Black Vinyl) Standard Anterior Trunk Support, Snap Closure	X-Large	24	12	14	8	2.5	8
VSX-V393	BioForm Non-Stretch (Black Vinyl) Custom Anterior Trunk Support, Snap Closure	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VSD-V393 becomes VSD-V322 when color Sun Yellow is desired.


If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.14.- BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support


Model	Description	Size	Dimension				
			A	B	C	D	E
VHS-V393	BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support	Small	10	14	6	8	2
VHM-V393	BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support	Medium	13	18	8	10	2.25
VHL-V393	BioForm Non-Stretch (Black Vinyl) “H”-Style Anterior Trunk Support	Large	16	21	10	12	2.5
VHX-V393	BioForm Non-Stretch (Black Vinyl) Custom “H”-Style Anterior Trunk Support	Custom	Template or Specs				

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VHD-V393 becomes VHD-V322 when color Sun Yellow is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.1.15.- BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support


Model	Description	Size	Dimension					
			A	B	C	D	E	F
VMT-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	XX-Small	10	5	5	4	1.75	1.25
VMD-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	X-Small	14	8	8.5	4	2	4
VMS-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	Small	17	9	10	5	2	5
VMM-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	Medium	20	10	11.5	6	2	6
VML-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	Large	22	11	12	7	2.5	7
VMO-V393	BioForm Non-Stretch (Black Vinyl) 6 Point Trunk Support	X-Large	24	12	14	8	2.5	8
VMX-V393	BioForm Non-Stretch (Black Vinyl) Custom 6 Point Trunk Support	Custom	Template or Specs					

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.
i.e.- VMD-V393 becomes VMD-V322 when color Sun Yellow is desired.

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-


“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.1.16.- BioForm Black Ballistic Nylon Standard Anterior Trunk Support




Model	Description	Size	Dimension					
			A	B	C	D	E	F
VNT	BioForm Ballistic Nylon Standard Anterior Trunk Support	XX-Small	10	5	5	4	1.75	1.25
VND	BioForm Ballistic Nylon Standard Anterior Trunk Support	X-Small	14	8	8.5	4	2	4
VNS	BioForm Ballistic Nylon Standard Anterior Trunk Support	Small	17	9	10	5	2	5
VNM	BioForm Ballistic Nylon Standard Anterior Trunk Support	Medium	20	10	11.5	6	2	6
VNL	BioForm Ballistic Nylon Standard Anterior Trunk Support	Large	22	11	12	7	2.5	7
VNO	BioForm Ballistic Nylon Standard Anterior Trunk Support	X-Large	24	12	14	8	2.5	8
VNX	BioForm Ballistic Nylon Custom Anterior Trunk Support	Custom	Template or Specs					

If Harness Guides are desired (to be mounted on a solid back insert), choose selection from section C 3.5-

B 2.2.- BioForm Abdominal Binders

The BioForm Abdominal Binders provide a crafty solution for proper positioning of a clients' mid section. It is constructed from a ½” piece of Sunmate foam covered in Neoprene with four straps attached to the body of the Abdominal Binder. Straps are equipped with Tri-Glides (for easy adjustment) and End fittings for easy mounting. Available in 3 sizes.

This option is a sensible choice when presented with the problem of a lordotic midsection that you would like to hold rearwards in the wheelchair.



BioForm Abdominal Binders

Model	Description	Size	Dimension	
			W	H
NAS-NB	BioForm Abdominal Binder, Four 1” Wide Straps, Fastex Buckle Closure, End fittings	Small	6	4
NAM-NB	BioForm Abdominal Binder, Four 1” Wide Straps, Fastex Buckle Closure, End fittings	Medium	9	5.5
NAL-NB	BioForm Abdominal Binder, Four 1” Wide Straps, Fastex Buckle Closure, End fittings	Large	12	7
NAX-NB	BioForm Custom Abdominal Binder	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- NAS-NB becomes NAS-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


B 2.3.- BioForm Positioning Belts

The BioForm positioning belts present a perfect complement to our seating systems. BioForm Positioning Belts are typically used at a 45 degree angle to properly position the pelvis or at a 90 degree angle across the femurs to maintain the leg contact with the seating surface.


The belts are available in 1", 1½" and 2" webbing widths and come with a choice of 3 closures; d-ring/velcro, fastex buckle and auto style buckle. Every belt comes standard with a tri-glide (for easy adjustment) and a d-ring /end fitting unit for easy mounting.

A padded belt can be created by choosing one of our padded sleeve options.


B 2.3.1.- BioForm Positioning Belts, Velcro Closure

	Model	Description	Dimension	
			W	L
	PBB	BioForm Positioning Belt, 1" webbing, D-Ring and Velcro Closure, End fittings	1	50
	PBC	BioForm Positioning Belt, 1 1/2" webbing, D-Ring and Velcro Closure, End fittings	1.5	60
	PBD	BioForm Positioning Belt, 2" webbing, D-Ring and Velcro Closure, End fittings	2	70
PBX	BioForm Custom Positioning Belt, webbing, D-Ring and Velcro Closure, End fittings	Template or Specs		

B 2.3.2.- BioForm Positioning Belts, Fastex Buckle Closure

	Model	Description	Dimension	
			W	L
	PEB	BioForm Positioning Belt, 1" webbing, Fastex Buckle Closure, End fittings	1	50
	PEC	BioForm Positioning Belt, 1 1/2" webbing, Fastex Buckle Closure, End fittings	1.5	60
	PED	BioForm Positioning Belt, 2" webbing, Fastex Buckle Closure, End fittings	2	70
PEX	BioForm Custom Positioning Belt, webbing, Fastex Buckle Closure, End fittings	Template or Specs		

B 2.3.3.- BioForm Positioning Belts, Auto Style Buckle Closure

	Model	Description	Dimension	
			W	L
	PFB	BioForm Positioning Belt, 1" webbing, Auto Style Buckle Closure, End fittings	1	50
	PFC	BioForm Positioning Belt, 1 1/2" webbing, Auto Style Buckle Closure, End fittings	1.5	60
	PFD	BioForm Positioning Belt, 2" webbing, Auto Style Buckle Closure, End fittings	2	70
PFX	BioForm Custom Positioning Belt, webbing, Auto Style Buckle Closure, End fittings	Template or Specs		

"If you don't see it... Contact Us"

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.4.- BioForm Positioning Belts with Split Pads


The BioForm Positioning Belts are also offered with “split pads” secured to both sides of the positioning belt. This option is effective for accurate protection of the A.S.I.S. bones as the pads can not slide off of the bony prominences.

The belts are available in 1”, 1½” and 2” webbing widths and come with a choice of 3 closures; d-ring/velcro, fastex buckle and auto style buckle. Every belt comes standard with a tri-glide (for easy adjustment) and a d-ring /end fitting unit for easy mounting.


B 2.4.1.- BioForm Positioning Belts, Split Pads, 2 each, Fastex Buckle Closure

	Model	Description	Dimension			
			Belt W	Belt L	Pad W	Pad L
	PIS	BioForm Positioning Belt, 1” wide webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	1	50	1.75	4
	PIT	BioForm Positioning Belt, 1½” wide webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	1.5	60	2.5	6
	PIU	BioForm Positioning Belt, 2” wide webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	2	70	4	8
	PIX	BioForm Custom Positioning Belt, webbing, Split Pads (2ea), Fastex Buckle Closure, End fittings	Template or Specs			

B 2.4.2.- BioForm Positioning Belts, Split Pads, 2 each, Auto-Style Buckle Closure

	Model	Description	Dimension			
			Belt W	Belt L	Pad W	Pad L
	PJV	BioForm Positioning Belt, 1” wide webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1	50	1.75	4
	PJW	BioForm Positioning Belt, 1½” wide webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1.5	60	2.5	6
	PJY	BioForm Positioning Belt, 2” wide webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	2	70	4	8
	PJX	BioForm Custom Positioning Belt, webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	Template or Specs			

B 2.4.3.- BioForm 4-Points Positioning Belts, Split Pads, 2 each, Auto-Style Buckle Closure

	Model	Description	Dimension			
			Belt W	Belt L	Pad W	Pad L
	PQV	BioForm 4-Points Positioning Belt, 1” wide webbing (¾” for anchoring), Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1	50	1.75	4
	PQW	BioForm 4-Points Positioning Belt, 1½” wide webbing (1” for anchoring), Split Pads (2ea), Auto-Style Buckle Closure, End fittings	1.5	60	2.5	6
	PQY	BioForm 4-Points Positioning Belt, 2” wide webbing (1” for anchoring), Split Pads (2ea), Auto-Style Buckle Closure, End fittings	2	70	4	8
	PQX	BioForm Custom 4-Points Positioning Belt, webbing, Split Pads (2ea), Auto-Style Buckle Closure, End fittings	Template or Specs			

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.5.- BioForm Chest Positioning Belts


The BioForm Chest Positioning Belts are beneficial for the consumer who needs assistance in keeping their upper trunk in contact to the solid back. These Positioning Belts can be installed to the Lateral Trunk Supports, to the Solid Back Inserts or to the wheelchair frame.

B 2.5.1.- BioForm Non-Stretch Chest Positioning Belts

The BioForm Non-Stretch Chest Positioning Belts are available in 1½” and 2” webbing widths. Both come with d-ring and velcro closure to allow for greater leverage when tightening. One side of the strap is placed in front of the client and passed through the d-ring on the other strap and brought back upon itself and secured via velcro. Positioning Belts have tri-glides for easy adjustment and end fittings for convenient mounting.

A padded Chest Positioning Belt can be created by adding one of the padded sleeve options.


BioForm Non-Stretch Chest Positioning Belts

	Model	Description	Dimension	
			W	L
	PKB	BioForm Chest Positioning Belt, 1½” wide webbing, D-Ring and Velcro Closure, End fittings	1.5	60
	PKC	BioForm Chest Positioning Belt, 2” wide webbing, D-Ring and Velcro Closure, End fittings	2	70
	PKX	BioForm Custom Chest Positioning Belt, webbing, D-Ring and Velcro Closure, End fittings	Template or Specs	

B 2.5.2.- BioForm Stretch (Black Neoprene) Chest Positioning Belts

The BioForm Stretch (Black Neoprene) Chest Positioning Belts are available with a pad of 8”, 11”, 14”, 17” and 20” in length. They come with velcro closure to allow for a quick setup. Each strap, short piece of webbing sewn to the Neoprene pads, overlap each other and attach via velcro in front of the client. Positioning Belts have tri-glides for easy adjustment and end fittings for convenient mounting. Available in x-small, small, medium, large and x-large. Custom sizes available by selecting the appropriate model #.

BioForm Stretch Chest Positioning Belts

	Model	Description	Size	Dimension (overlap)	
				W	L
	PDD-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 1.5” End fittings	X-Small	2.5	6
	PDS-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	Small	4	9
	PDM-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	Medium	4	12
	PDL-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	Large	4	15
	PDO-NB	BioForm Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, 2” End fittings	X-Large	5	20
	PDX-NB	BioForm Custom Stretch (Black Neoprene) Chest Positioning Belt, Velcro Closure, End fittings	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- PDD-NB becomes PDD-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.6.- BioForm Calf Straps

The BioForm Calf Strap is a simple solution to assist in holding the lower extremities from slipping rearwards between the legrests. The strap is 2” wide webbing with a velcro closure. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.

The Calf Strap can be made padded by adding one of the BioForm Padded Sleeve options.

BioForm Calf Straps


	Model	Description	Size	Dimension	
				W	L
	PCS	BioForm Calf Strap, 2” wide webbing with a Velcro closure	Small		30
	PCM	BioForm Calf Strap, 2” wide webbing with a Velcro closure	Medium		35
	PCL	BioForm Calf Strap, 2” wide webbing with a Velcro closure	Large		40
	PCX	BioForm Custom Calf Strap, webbing with a Velcro closure	Custom		Template or Specs

B 2.7.- BioForm Thigh Straps

The BioForm Thigh Straps are typically chosen when a belt positioned at 90 degrees to the thighs is useful in keeping the femurs in contact with the seat cushion. The straps are available in 1½” and 2” wide webbing with D-Ring and Velcro Closure.

The Thigh Strap can be made padded by adding one of the BioForm Padded Sleeve options.

BioForm Thigh Straps

	Model	Description	Dimension	
			W	L
	PTB	BioForm Thigh Strap, 1½” wide webbing, D-Ring and Velcro Closure, End fittings	1.5	60
	PTC	BioForm Thigh Strap, 2” wide webbing, D-Ring and Velcro Closure, End fittings	2	70
	PTX	BioForm Custom Thigh Strap, webbing, D-Ring and Velcro Closure, End fittings		Template or Specs

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

B 2.8.- BioForm Padded Sleeves



The BioForm Padded Sleeves are the perfect selection when any belt currently in use produces excessive pressure areas. They can be retrofitted to any positioning belt.

The BioForm line of Padded Sleeves are made from ½” of Sunmate foam covered in Neoprene (Rubatex). Two attaching loops are sewn onto pad. One side is fixed to allow belt or strap to pass through and the other is a Velcro overlap to allow for quick and easy removal so belt or strap can be easily undone. This securing method assures caregivers that the pad will stay attached to belt or strap while still easily removable for transfers. Pads are available to fit our 1 ½” belts and our 2” belts and come in 10”, 12” and 14” lengths. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.

B 2.8.1.- BioForm Black Padded Sleeves For 2” Belts

Model	Description	Size	Dimension	
			L	W
PGS-NB	BioForm Black Padded Sleeve for 2” Belt	Small	10”	
PGM-NB	BioForm Black Padded Sleeve for 2” Belt	Medium	12.5”	
PGL-NB	BioForm Black Padded Sleeve for 2” Belt	Large	15”	
PGX-NB	BioForm Custom Black Padded Sleeve for 2” Belt	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- PGM-NB becomes PGM-NE when color Eggplant is desired.

B 2.8.2.- BioForm Black Padded Sleeves For 1½” Belts

Model	Description	Size	Dimension	
			L	W
PHS-NB	BioForm Black Padded Sleeve for 1½” Belt	Small	10”	
PHM-NB	BioForm Black Padded Sleeve for 1½” Belt	Medium	12.5”	
PHL-NB	BioForm Black Padded Sleeve for 1½” Belt	Large	15”	
PHX-NB	BioForm Custom Black Padded Sleeve for 1½” Belt	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section. i.e.- PHM-NB becomes PHM-NE when color Eggplant is desired.

“If you don’t see it... Contact Us”

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – www.BioDynamics.us

E-Mail: quote@BioDynamics.us Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.