

## B 3.- Lower Extremity Supports

The Lower Extremity Supports offered in this section enable proper positioning of the lower extremities. Every selection comes with a mounting hardware option.


### B 3.1.- BioForm Foot Positioners



The BioForm line of foot and ankle positioners include: Ankle Positioners, Footplate Extenders, Ankle Positioners with adjustable Plate and Foot Positioners,

- **BioForm Ankle Positioners** are a 3 point ankle cuff forming a cup around the ankle to maintain proper positioning of the foot/ankle on a footplate. They can be mounted to any slotted footplate or our BioForm Multi-Slots Footplate Extenders. All mounting hardware is included. Fastex Buckle closure. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.
- **BioForm Multi-Slots Footplate Extenders** allow for the foot to be positioned in a location off of the footplate; anterior, posterior or angled. To these multi-slots extenders, our foot and ankle positioning accessories can be mounted in the multiple mounting slots available. All mounting hardware is included. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.
- **BioForm Ankle Positioners with adjustable Plate** includes a BioForm Footplate Extender and a BioForm Ankle Positioner allowing for a proper positioning of the foot/ankle in a location of the footplate. All mounting hardware is included. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.
- **BioForm Foot Positioners** provide for maximum positioning capability of the foot and ankle. They consist of a BioForm Multi-Slots Footplate Extender, BioForm Ankle Positioner and a BioForm Padded Ankle/Toe Strap. All mounting hardware is included. Sold in pairs. Available in small, medium and large. Custom sizes available by selecting the appropriate model.

#### B 3.1.1.- BioForm Ankle Positioners

	Model	Description	Size	Fits Ankle Circumference	
				From	To
	FAS-NB	BioForm Ankle Positioner, pair	Small	5	9
	FAM-NB	BioForm Ankle Positioner, pair	Medium	7.5	12
	FAL-NB	BioForm Ankle Positioner, pair	Large	10	15
	FAX-NB	BioForm Custom Ankle Positioner	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- FAS-NB becomes FAS-NE when color Eggplant is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.1.2.- BioForm Multi-Slot Footplate Extenders




Model	Description	Size	Dimension	
			W	L
FXS	BioForm Multi-Slot Footplate Extenders, pair	Small	4	8
FXM	BioForm Multi-Slot Footplate Extenders, pair	Medium	5	9.5
FXL	BioForm Multi-Slot Footplate Extenders, pair	Large	5	11
FXX	BioForm Custom Multi-Slot Footplate Extenders	Custom	Template or Specs	

### B 3.1.3.- BioForm Ankle Positioner with Adjustable Plate

	Model	Description	Size	Fits Ankle Circumference		Dimension (Plate)	
				From	To	W	L
	FBS-NB	BioForm Ankle Positioner With Adjustable Plate, pair	Small	5	9	4	8
	FBM-NB	BioForm Ankle Positioner With Adjustable Plate, pair	Medium	7.5	12	5	9.5
	FBL-NB	BioForm Ankle Positioner With Adjustable Plate, pair	Large	10	15	5	11
	FBX-NB	BioForm Custom Ankle Positioner With Adjustable Plate	Custom	Template or Specs			

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- FBS-NB becomes FBS-NE when color Eggplant is desired.

### B 3.1.4.- BioForm Foot Positioner

	Model	Description	Size	Fits Ankle Circumference		Dimension (Plate)	
				From	To	W	L
	FPS-NB	BioForm Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender) , pair	Small	5	9	4	8
	FPM-NB	BioForm Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender) , pair	Medium	7.5	12	5	9.5
	FPL-NB	BioForm Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender) , pair	Large	10	15	5	11
	FPX-NB	BioForm Custom Foot Positioner (Ankle Positioner, Toe Strap and Footplate Extender)	Custom	Template or Specs			

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- FPS-NB becomes FPS-NE when color Eggplant is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

## B 3.2.- BioForm Footplates


### B 3.2.1.- BioForm One Piece, Footplates

BioDynamics manufactures One Piece, Footplates that are shipped with pre-drilled holes (countersunk) and mounting hardware for easy mounting in the appropriate position. Our One Piece, Footplates are made with a ½” High Density Polyethylene (H.D.P.E.) base. The padded version has a ¼” of H.D.P.E., summate foam and Vinyl (Naugahyde) covered module secured to the base with hardware.

Padded and non-padded One Piece, Footplates are available with fixed hardware or flip up hardware.. The flip up hardware allows the footrests to be swung away/removed when necessary. When the footrests are put back on the chair, the One Piece, Footplate swings down to the opposite footplate from where it is attached and has an indexing pin that will set into a hole the RTS will drill into that footplate. This makes the One Piece, Footplate more secure when in operation.


Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#.

#### B 3.2.1.1.- BioForm One Piece, Fixed Footplate



Model	Description	Dimension	
		W	D
CFS	BioForm One Piece, Fixed Footplate, for 14” Wheelchairs	11.5	8
CFM	BioForm One Piece, Fixed Footplate, for 16” Wheelchairs	13.5	9
CFL	BioForm One Piece, Fixed Footplate, for 18” Wheelchairs	15.5	10
CFO	BioForm One Piece, Fixed Footplate, for 20” Wheelchairs	17.5	10
CFX	BioForm Custom One Piece, Fixed Footplate	Template or Specs	

#### B 3.2.1.2.- BioForm Padded One Piece, Fixed Footplate



Model	Description	Dimension	
		W	D
CGS-V393	BioForm Padded One Piece, Fixed Footplate, for 14” Wheelchairs	11.5	8
CGM-V393	BioForm Padded One Piece, Fixed Footplate, for 16” Wheelchairs	13.5	9
CGL-V393	BioForm Padded One Piece, Fixed Footplate, for 18” Wheelchairs	15.5	10
CGO-V393	BioForm Padded One Piece, Fixed Footplate, for 20” Wheelchairs	17.5	10
CGX-V393	BioForm Padded Custom One Piece, Fixed Footplate	Template or Specs	

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- CGS-V393 becomes CGS-V354 when Vinyl Royal Blue is desired.


***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.2.1.3.- BioForm One Piece, Flip-Up Footplate



Model	Description	Dimension	
		W	D
CES	BioForm One Piece, Flip-Up Footplate, Self Centering Locking, for 14" Wheelchairs	11.5	8
CEM	BioForm One Piece, Flip-Up Footplate, Self Centering Locking,, for 16" Wheelchairs	13.5	9
CEL	BioForm One Piece, Flip-Up Footplate, Self Centering Locking, for 18" Wheelchairs	15.5	10
CEO	BioForm One Piece, Flip-Up Footplate, Self Centering Locking, for 20" Wheelchairs	17.5	10
CEX	BioForm Custom One Piece, Flip-Up Footplate, Self Centering Locking,	Template or Specs	

### B 3.2.1.4.- BioForm Padded One Piece, Flip-Up Footplate



Model	Description	Dimension	
		W	D
CDS-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 14" Wheelchairs	11.5	8
CDM-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 16" Wheelchairs	13.5	9
CDL-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 18" Wheelchairs	15.5	10
CDO-V393	BioForm Padded One Piece, Flip-Up Footplate, Self Centering Locking, for 20" Wheelchairs	17.5	10
CDX-V393	BioForm Custom Padded One Piece, Flip-Up Footplate, Self Centering Locking,	Template or Specs	

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- CDS-V393 becomes CDS-V354 when Vinyl Royal Blue is desired.

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.2.2.- BioForm Footplates Extenders

BioDynamics manufactures Footplates Extenders that are shipped with pre-drilled holes (countersunk) and mounting hardware for easy mounting in the appropriate position. Our Footplates Extenders are made with a ½” High Density Polyethylene (H.D.P.E.) base. The padded version has a ¼” of H.D.P.E., summate foam and Vinyl (Naugahyde) covered module secured to the base with hardware.


Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#.

#### B 3.2.2.1.- BioForm Footplates Extender



Model	Description	Dimension	
		W	D
CIS	BioForm Footplates Extender, for 14” Wheelchairs, pair.	5.5	9
CIM	BioForm Footplates Extender, for 16” Wheelchairs, pair.	6.5	10
CIL	BioForm Footplates Extender, for 18” Wheelchairs, pair.	7.5	11
CIO	BioForm Footplates Extender, for 20” Wheelchairs, pair.	8.5	12
CIX	BioForm Custom Footplates Extender	Template or Specs	

#### B 3.2.2.2.- BioForm Padded Footplates Extender



Model	Description	Dimension	
		W	D
CJS-V393	BioForm Padded Footplates Extender, for 14” Wheelchairs, pair.	5.5	9
CJM-V393	BioForm Padded Footplates Extender, for 16” Wheelchairs, pair.	6.5	10
CJL-V393	BioForm Padded Footplates Extender, for 18” Wheelchairs, pair.	7.5	11
CJO-V393	BioForm Padded Footplates Extender, for 20” Wheelchairs, pair.	8.5	12
CJX-V393	BioForm Custom Padded Footplates Extender	Template or Specs	

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- CJS-V393 becomes CJS-V354 when Vinyl Royal Blue is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.2.3.- BioForm Footplate Pads

#### B 3.2.3.1.- BioForm Soft Footplate Pads

BioDynamics manufactures BioForm Soft Footplate Pads that are designed to protect the consumers feet from the footplates. The Footplate Pads have a pocket sewn onto the underside of the pad which enables the entire pad to slip onto the footplate. The pad is then further secured to the footplate with a velcro strap, also, sewn onto the pad. The pads are made of ½” of Sunmate foam covered in Neoprene.

Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#. Sold in pairs.



**BioForm Soft Footplate Pads**

Model	Description	Dimension	
		W	D
CMS-NB	BioForm Soft Footplate Pads, for 14” Wheelchairs, pair	11.5	8
CMM-NB	BioForm Soft Footplate Pads, for 16” Wheelchairs, pair	13.5	9
CML-NB	BioForm Soft Footplate Pads, for 18” Wheelchairs, pair	15.5	10
CMO-NB	BioForm Soft Footplate Pads, for 20” Wheelchairs, pair	17.5	10
CMX-NB	BioForm Custom Soft Footplate Pads	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- CMS-NB becomes CMS-NE when color Eggplant is desired.

#### B 3.2.3.2.- BioForm Footplate Pads

BioDynamics manufactures Footplate Pads that are designed to protect the consumers feet from the footplates. These Footplate Pads are secured to the footplate by inserting a pair of bolts through the underside of the footplate and into the T-nuts embedded into the plastic that is on the underside of the Footplate Pad. The pads are made of ½” of Sunmate foam covered in Neoprene which is secured to a ¼” plastic (HDPE) base.

Our standard sizes are available to fit 14”, 16”, 18” and 20” wheelchairs. If a custom size is needed, select the appropriate model#. Sold in pairs.



**BioForm Footplate Pads**

Model	Description	Size	Dimension	
			W	D
CNS-NB	BioForm Footplate Pads, pair	Small	4	6.5
CNM-NB	BioForm Footplate Pads, pair	Medium	5	6.5
CNL-NB	BioForm Footplate Pads, pair	Large	6	8
CNO-NB	BioForm Footplate Pads, pair	X-Large	7	8
CNX-NB	BioForm Custom Footplate Pads	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- CNS-NB becomes CNS-NE when color Eggplant is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



### B 3.3.- BioForm Sandals and Straps



Our BioForm Sandals securely hold the feet in proper position utilizing either: an adjustable Fastex Buckle strap, adjustable Velcro strap or an adjustable leather strap. The bottom of the sandal is slotted for ease in finding the proper mounting location. Sandals are available with either one or two straps and are sold in pairs. All mounting hardware is included. Available in small, medium and large. Custom sizes available by selecting the appropriate model #.

#### B 3.3.1.- BioForm Sandals with one Strap w/Velcro



Model	Description	Size	Dimension	
			W	D
LIS	BioForm Sandals with one Strap w/Velcro, pair	Small	7	4
LIM	BioForm Sandals with one Strap w/Velcro, pair	Medium	9	4.5
LIL	BioForm Sandals with one Strap w/Velcro, pair	Large	11	5
LIX	BioForm Custom Sandals with one Strap w/Velcro	Custom	Template or Specs	

#### B 3.3.2.- BioForm Sandals with one Strap w/Buckle



Model	Description	Size	Dimension	
			W	D
LJS	BioForm Sandals with one Strap w/Buckle, pair	Small	7	4
LJM	BioForm Sandals with one Strap w/Buckle, pair	Medium	9	4.5
LJL	BioForm Sandals with one Strap w/Buckle, pair	Large	11	5
LJX	BioForm Custom Sandals with one Strap w/Buckle	Custom	Template or Specs	


***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.3.3.- BioForm Sandals with one Leather Strap




Price	Model	Description	Size	Dimension	
				W	D
225	LMS	BioForm Sandals with one Leather Strap, pair	Small	7	4
225	LMM	BioForm Sandals with one Leather Strap, pair	Medium	9	4.5
225	LML	BioForm Sandals with one Leather Strap, pair	Large	11	5
Call	LMX	BioForm Custom Sandals with one Leather Strap	Custom	Template or Specs	

### B 3.3.4.- BioForm Sandals with two Straps w/Velcro



Model	Description	Size	Dimension	
			W	D
LKS	BioForm Sandals with two Straps w/Velcro, pair	Small	7	4
LKM	BioForm Sandals with two Straps w/Velcro, pair	Medium	9	4.5
LKL	BioForm Sandals with two Straps w/Velcro, pair	Large	11	5
LKX	BioForm Custom Sandals with two Straps w/Velcro	Custom	Template or Specs	

### B 3.3.5.- BioForm Sandals with two Straps w/Buckle



Model	Description	Size	Dimension	
			W	D
LLS	BioForm Sandals with two Straps w/Buckle, pair	Small	7	4
LLM	BioForm Sandals with two Straps w/Buckle, pair	Medium	9	4.5
LLL	BioForm Sandals with two Straps w/Buckle, pair	Large	11	5
LLX	BioForm Custom Sandals with two Straps w/Buckle	Custom	Template or Specs	

***“If you don’t see it.... Contact Us”***


BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



### B 3.3.6.- BioForm Sandals with two Leather Straps



Model	Description	Size	Dimension	
			W	D
LNS	BioForm Sandals with two Leather Straps, pair	Small	7	4
LNМ	BioForm Sandals with two Leather Straps, pair	Medium	9	4.5
LNL	BioForm Sandals with two Leather Straps, pair	Large	11	5
LNX	BioForm Custom Sandals with two Leather Straps	Custom	Template or Specs	


### B 3.4.- Solid Heel Cups

#### B 3.4.1.- BioForm Solid Heel Cups with Straps w/ Velcro

This choice for foot positioning is optimal when the primary goal is to secure the foot in a specific location on the footplate while allowing for some front of foot movement. Heel Cups are very beneficial for consumers with significant lower extremity tone.

BioForm Solid Heel Cups function in a similar fashion as the BioForm Sandals do. The base of the Heelcup is made from ½" HDPE, pre-drilled with countersunk holes for seamless mounting. The rear wall is made from ¼" kydex and is attached to the Heel Cup base with screws. The strap is affixed to the base of the Heelcup and is available with a Fastex Buckle closure or a Velcro closure.

Available in small, medium and large and come complete with all mounting hardware. Custom sizes available by selecting the appropriate model #. Sold in pairs.



**BioForm Solid Heel Cups with Straps w/ Velcro**

Model	Description	Size	Dimension		
			W	D	H
LAS	BioForm Solid Heel Cups with Straps w/ Velcro, pair	Small	4	7	2.5
LAM	BioForm Solid Heel Cups with Straps w/ Velcro, pair	Medium	4.5	9	2.5
LAL	BioForm Solid Heel Cups with Straps w/ Velcro, pair	Large	5	11	2.5
LAX	BioForm Custom Solid Heel Cups with Straps w/ Velcro	Custom	Template or Specs		


***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.4.2.- BioForm Solid Heel Cups with Straps w/ Fastex Buckle




Model	Description	Size	Dimension		
			W	D	H
LBS	BioForm Solid Heel Cups with Straps w/ Fastex Buckle, pair	Small	4	7	2.5
LBM	BioForm Solid Heel Cups with Straps w/ Fastex Buckle, pair	Medium	4.5	9	2.5
LBL	BioForm Solid Heel Cups with Straps w/ Fastex Buckle, pair	Large	5	11	2.5
LBX	BioForm Custom Solid Heel Cups with Straps w/ Fastex Buckle	Custom	Template or Specs		

### B 3.5.- BioForm Foot Support Accessories

#### B 3.5.1.- BioForm Scuba Heel Loops



BioForm Scuba Heel Loops are an industry first. These “crush proof” heel loops are the answer to the chronic issue of heel loops getting crushed to the surface of a footplate. They are fabricated from a special webbing and come complete with all of the necessary mounting hardware. Available in small, medium and large and are sold in pairs. Custom sizes available by selecting the appropriate model #.



**BioForm Scuba Heel Loops**

Model	Description	Size	Dimension	
			L	H
LCS	BioForm Scuba Heel Loops, pair	Small	6	2
LCM	BioForm Scuba Heel Loops, pair	Medium	8	2
LCL	BioForm Scuba Heel Loops, pair	Large	10	2
LCX	BioForm Custom Scuba Heel Loops	Custom	Template or Specs	

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)


E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


### B 3.5.2.- BioForm Ankle or Toe Straps

The BioForm Ankle or Toe Straps are an effective means of holding a foot or ankle in place upon a footplate. They are available in 2 styles: padded and non-padded. Also, the non-padded style comes with the choice of Velcro or Fastex Buckle closure. Each strap is 1" wide, has Tri-Glides for easy adjusting and comes with D-Rings that house the end fittings. Our end fittings provide a versatile tool for easy and convenient mounting. Available in small, medium and large and are sold in pairs. Custom sizes available by selecting the appropriate model #.


#### B 3.5.2.1.- BioForm Ankle or Toe Straps, Velcro Closure

	Model	Description	Size	Dimension	
				From	To
	MDS	BioForm Ankle or Toe Straps, 1", D-Ring Velcro Closure, End fittings, pair	Small	7	10.5
	MDM	BioForm Ankle or Toe Straps, 1", D-Ring Velcro Closure, End fittings, pair	Medium	8	13.5
	MDL	BioForm Ankle or Toe Straps, 1", D-Ring Velcro Closure, End fittings, pair	Large	9	18.5
	MDX	BioForm Custom Ankle or Toe Straps, D-Ring Velcro Closure, End fittings	Custom	Template or Specs	

#### B 3.5.2.2.- BioForm Ankle or Toe Straps, Fastex Buckle Closure

	Model	Description	Size	Dimension	
				From	To
	MFS	BioForm Ankle or Toe Straps, 1", Fastex Buckle Closure, End fittings, pair	Small	7	10.5
	MFM	BioForm Ankle or Toe Straps, 1", Fastex Buckle Closure, End fittings, pair	Medium	8	13.5
	MFL	BioForm Ankle or Toe Straps, 1", Fastex Buckle Closure, End fittings, pair	Large	9	18.5
	MFY	BioForm Custom Ankle or Toe Straps, Fastex Buckle Closure, End fittings	Custom	Template or Specs	

#### B 3.5.2.3.- BioForm Padded Ankle or Toe Straps, Fastex Buckle Closure

	Model	Description	Size	Dimension	
				From	To
	MPS-NB	BioForm Padded Ankle or Toe Straps, 1", Fastex Buckle Closure, End fittings, pair	Small	7	10.5
	MPM-NB	BioForm Padded Ankle or Toe Straps, 1", Fastex Buckle Closure, End fittings, pair	Medium	8	13.5
	MPL-NB	BioForm Padded Ankle or Toe Straps, 1", Fastex Buckle Closure, End fittings, pair	Large	9	18.5
	MPX-NB	BioForm Custom Padded Ankle or Toe Straps, Fastex Buckle Closure, End fittings	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- MPS-NB becomes MPS-NE when color Eggplant is desired.

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)


E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

#### B 3.5.2.4.- Pads for Ankle or Toe Straps

The BioForm Pads for Ankle or Toe Straps provide an opportunity to add padded protection for straps that were ordered without padding. The pads are sleeve type for easy installation.

*BioForm Padded Sleeves for Ankle or Toe Straps*


	Model	Description	Dimension	
			W	L
	LSP-NB	BioForm Padded Sleeves for Ankle or Toe Straps, pair	2	5.5
	LSX-NB	BioForm Custom Padded Sleeves for Ankle or Toe Straps	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- LSP-NB becomes LSP-NE when color Eggplant is desired.

#### B 3.5.2.5.- BioForm Leather Straps

The BioForm Leather Straps are a versatile tool for very aggressive positioning needs of the foot, ankle or any other area with strong, unwanted movements. They are sold in pairs and are 1" in width. The straps are 2 pieces that are secured together with a Standard Duty Steel "belt type" buckle and have a hole at both ends for easy mounting.

**BioForm Leather Straps**

	Model	Description	Size	Dimension	
				From	To
	MMS	BioForm Leather Ankle or Toe Straps, 1" wide, pair	Small	12	16
	MMM	BioForm Leather Ankle or Toe Straps, 1" wide, pair	Medium	14	20
	MML	BioForm Leather Ankle or Toe Straps, 1" wide, pair	Large	19	24


***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.6.- BioForm Footboxes





The BioForm line of Footboxes offer a wide array of styles and mounting options. Our 2 styles of Footboxes are: single leg boxes (designed to accommodate one leg per Footbox), and a 1 piece Footbox (designed to accommodate both legs in one Footbox). The Footboxes are available with hardware for: mounting to footplates, mounting to hanger tubes or mounting to a solid seat pan. We also offer our Footboxes without mounting hardware to offer you a more local/custom mounting solution. All Footbox options are available padded or non-padded. Padded units have the pad modules secured to the walls of the Footboxes with t-nuts and screws for added durability.

Our standard sizes are available to fit 14", 16", 18" and 20" wheelchairs. If a custom size is needed, select the appropriate model#.

#### B 3.6.1.- BioForm, Double Leg (1 piece), Non-padded Footboxes

##### B 3.6.1.1.- BioForm, Double Leg (1 piece), Non-padded Footboxes, No Hardware



Model	Description	Dimension		
		W	D	H
EAS	BioForm Double Leg Footbox, No Hardware, for 14" Wheelchairs	11.5	8	8
EAM	BioForm Double Leg Footbox, No Hardware, for 16" Wheelchairs	13.5	9	9
EAL	BioForm Double Leg Footbox, No Hardware, for 18" Wheelchairs	15.5	10	10.5
EAO	BioForm Double Leg Footbox, No Hardware, for 20" Wheelchairs	17.5	11	12
EAX	BioForm Custom Double Leg Footbox, No Hardware	Template or Specs		



***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345



© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

**B 3.6.1.2.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Nuts and Bolts to secure from footplates.**



Model	Description	Dimension		
		W	D	H
EBS	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs	11.5	8	8
EBM	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs	13.5	9	9
EBL	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs	15.5	10	10.5
EBO	BioForm Double Leg Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs	17.5	11	12
EBX	BioForm Custom Double Leg Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

**B 3.6.1.3.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Hardware to secure from hangers.**



Model	Description	Dimension		
		W	D	H
ECS	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 14" Wheelchairs	11.5	8	8
ECM	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 16" Wheelchairs	13.5	9	9
ECL	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 18" Wheelchairs	15.5	10	10.5
ECO	BioForm Double Leg Footbox, Hardware to secure from Hangers, for 20" Wheelchairs	17.5	11	12
ECX	BioForm Custom Double Leg Footbox, Hardware to secure from Hangers	Template or Specs		

***"If you don't see it.... Contact Us"***



BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.




**B 3.6.1.4.- BioForm, Double Leg (1 piece), Non-padded Footboxes, Hardware to secure from seat pan.**



Model	Description	Dimension		
		W	D	H
EDS	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs	11.5	8	8
EDM	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs	13.5	9	9
EDL	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs	15.5	10	10.5
EDO	BioForm Double Leg Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs	17.5	11	12
EDX	BioForm Custom Double Leg Footbox, Hardware to secure from Seat Pan	Template or Specs		

**B 3.6.2.- BioForm, Single Leg (2 pieces), Non-padded Footboxes**

**B 3.6.2.1.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, No Hardware**



Model	Description	Dimension		
		W	D	H
EES	BioForm Single Leg Footbox, No Hardware, for 14" Wheelchairs (pair)	5.5	8	8
EEM	BioForm Single Leg Footbox, No Hardware, for 16" Wheelchairs (pair)	6.5	9	9
EEL	BioForm Single Leg Footbox, No Hardware, for 18" Wheelchairs (pair)	7.5	10	10.5
EEO	BioForm Single Leg Footbox, No Hardware, for 20" Wheelchairs (pair)	8.5	11	12
EEX	BioForm Custom Single Leg Footbox, No Hardware	Template or Specs		


***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345



© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

**B 3.6.2.2.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Nuts and Bolts to secure from footplates.**



Model	Description	Dimension		
		W	D	H
EFS	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs (pair)	5.5	8	8
EFM	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs (pair)	6.5	9	9
EFL	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs (pair)	7.5	10	10.5
EFO	BioForm Single Leg Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs (pair)	8.5	11	12
EFX	BioForm Custom Single Leg Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

**B 3.6.2.3.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Hardware to secure from hangers.**



Model	Description	Dimension		
		W	D	H
EGS	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 14" Wheelchairs (pair)	5.5	8	8
EGM	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 16" Wheelchairs (pair)	6.5	9	9
EGL	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 18" Wheelchairs (pair)	7.5	10	10.5
EGO	BioForm Single Leg Footbox, Hardware to secure from Hangers, for 20" Wheelchairs (pair)	8.5	11	12
EGX	BioForm Custom Single Leg Footbox, Hardware to secure from Hangers	Template or Specs		

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

**B 3.6.2.4.- BioForm, Single Leg (2 pieces), Non-padded Footboxes, Hardware to secure from seat pan.**



Model	Description	Dimension		
		W	D	H
EHS	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs (pair)	5.5	8	8
EHM	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs (pair)	6.5	9	9
EHL	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs (pair)	7.5	10	10.5
EHO	BioForm Single Leg Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs (pair)	8.5	11	12
EHX	BioForm Custom Single Leg Footbox, Hardware to secure from Seat Pan	Template or Specs		

**B 3.6.3.- BioForm, Double Leg (1 piece), Padded Footboxes**

**B 3.6.3.1.- BioForm, Double Leg (1 piece), Padded Footboxes, No Hardware**



Model	Description	Dimension		
		W	D	H
EIS-V393	BioForm Double Leg Padded Footbox, No Hardware, for 14" Wheelchairs	11.5	8	8
EIM-V393	BioForm Double Leg Padded Footbox, No Hardware, for 16" Wheelchairs	13.5	9	9
EIL-V393	BioForm Double Leg Padded Footbox, No Hardware, for 18" Wheelchairs	15.5	10	10.5
EIO-V393	BioForm Double Leg Padded Footbox, No Hardware, for 20" Wheelchairs	17.5	11	12
EIX-V393	BioForm Custom Double Leg Padded Footbox, No Hardware	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- EIS-V393 becomes EIS-V354 when Vinyl Royal Blue is desired.


***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.


**B 3.6.3.2.- BioForm, Double Leg (1 piece), Padded Footboxes, Nuts and Bolts to secure from footplates.**



Model	Description	Dimension		
		W	D	H
EJS-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs	11.5	8	8
EJM-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs	13.5	9	9
EJL-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs	15.5	10	10.5
EJO-V393	BioForm Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs	17.5	11	12
EJX-V393	BioForm Custom Double Leg Padded Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- EJS-V393 becomes EJS-V354 when Vinyl Royal Blue is desired.

**B 3.6.3.3.- BioForm, Double Leg (1 piece), Padded Footboxes, Hardware to secure from hangers.**



Model	Description	Dimension		
		W	D	H
EKS-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 14" Wheelchairs	11.5	8	8
EKM-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 16" Wheelchairs	13.5	9	9
EKL-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 18" Wheelchairs	15.5	10	10.5
EKO-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Hangers, for 20" Wheelchairs	17.5	11	12
EKX-V393	BioForm Custom Double Leg Padded Footbox, Hardware to secure from Hangers	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- EKS-V393 becomes EKS-V354 when Vinyl Royal Blue is desired.



***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

#### B 3.6.3.4.- BioForm, Double Leg (1 piece), Padded Footboxes, Hardware to secure from seat pan.





Model	Description	Dimension		
		W	D	H
ELS-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs	11.5	8	8
ELM-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs	13.5	9	9
ELL-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs	15.5	10	10.5
ELO-V393	BioForm Double Leg Padded Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs	17.5	11	12
ELX-V393	BioForm Custom Double Leg Padded Footbox, Hardware to secure from Seat Pan	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- ELS-V393 becomes ELS-V354 when Vinyl Royal Blue is desired.

#### B 3.6.4.- BioForm, Single Leg (2 pieces), Padded Footboxes

##### B 3.6.4.1.- BioForm, Single Leg (2 pieces), Padded Footboxes, No Hardware



Model	Description	Dimension		
		W	D	H
EMS-V393	BioForm Single Leg Padded Footbox, No Hardware, for 14" Wheelchairs (pair)	5.5	8	8
EMM-V393	BioForm Single Leg Padded Footbox, No Hardware, for 16" Wheelchairs (pair)	6.5	9	9
EML-V393	BioForm Single Leg Padded Footbox, No Hardware, for 18" Wheelchairs (pair)	7.5	10	10.5
EMO-V393	BioForm Single Leg Padded Footbox, No Hardware, for 20" Wheelchairs (pair)	8.5	11	12
EMX-V393	BioForm Custom Single Leg Padded Footbox, No Hardware	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- EMS-V393 becomes EMS-V354 when Vinyl Royal Blue is desired.

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

**B 3.6.4.2.- BioForm, Single Leg (2 pieces), Padded Footboxes, Nuts and Bolts to secure from footplates.**



Model	Description	Dimension		
		W	D	H
ENS-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 14" Wheelchairs (pair)	5.5	8	8
ENM-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 16" Wheelchairs (pair)	6.5	9	9
ENL-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 18" Wheelchairs (pair)	7.5	10	10.5
ENO-V393	BioForm Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates, for 20" Wheelchairs (pair)	8.5	11	12
ENX-V393	BioForm Custom Single Leg Padded Footbox, Nuts and Bolts to secure from Footplates	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- ENS-V393 becomes ENS-V354 when Vinyl Royal Blue is desired.

**B 3.6.4.3.- BioForm, Single Leg (2 pieces), Padded Footboxes, Hardware to secure from hangers.**



Model	Description	Dimension		
		W	D	H
EOS-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 14" Wheelchairs (pair)	5.5	8	8
EOM-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 16" Wheelchairs (pair)	6.5	9	9
EOL-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 18" Wheelchairs (pair)	7.5	10	10.5
EOO-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Hangers, for 20" Wheelchairs (pair)	8.5	11	12
EOX-V393	BioForm Custom Single Leg Padded Footbox, Hardware to secure from Hangers	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- EOS-V393 becomes EOS-V354 when Vinyl Royal Blue is desired.

***"If you don't see it.... Contact Us"***



BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.



**B 3.6.4.4.- BioForm, Single Leg (2 pieces), Padded Footboxes, Hardware to secure from seat pan.**



Model	Description	Dimension		
		W	D	H
EPS-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 14" Wheelchairs (pair)	5.5	8	8
EPM-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 16" Wheelchairs (pair)	6.5	9	9
EPL-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 18" Wheelchairs (pair)	7.5	10	10.5
EPO-V393	BioForm Single Leg Padded Footbox, Hardware to secure from Seat Pan, for 20" Wheelchairs (pair)	8.5	11	12
EPX-V393	BioForm Custom Single Leg Padded Footbox, Hardware to secure from Seat Pan	Template or Specs		

If covering other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section.  
i.e.- EPS-V393 becomes EPS-V354 when Vinyl Royal Blue is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

## B 3.7.- BioForm Calf Support

### B 3.7.1.- BioForm Soft Calf Panels

BioForm Soft Calf Panels are a useful option when assistance is needed in keeping the consumers' legs from sliding rearwards off of the footplates.

The line of BioForm Soft Calf Panels are available in 2 styles: Non-Stretch, Vinyl (Nauahyde) for a more firm support and Stretch, double layer of Neoprene for a panel that is softer and has more flexibility. Each comes with straps for easy securing to the footrest/legrest assembly and close with Velcro or a Fastex Buckle.

Available in small, medium, large and x-large. Custom sizes available by selecting the appropriate model #.

Each style is available in 4 sizes.

#### B 3.7.1.1.- BioForm Stretch Soft Black Calf Panels with Fastex Buckle Closure



Model	Description	Size	Dimension	
			W	H
NCS-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	Small	16	7
NCM-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	Medium	18	8
NCL-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	Large	20	8
NCO-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 2 Straps w/ Fastex Buckles	X-Large	22	9
NCX-NB	BioForm Custom Stretch Soft Black Calf Panel, Layers of Neoprene, Fastex Buckles	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- NCS-NB becomes NCS-NE when color Eggplant is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.7.1.2.- BioForm Stretch Soft Black Calf Panels with Velcro Closure (4+2 straps)



Model	Description	Size	Dimension	
			W	H
NDS-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	Small	16	7
NDM-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	Medium	18	8
NDL-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	Large	20	8
NDO-NB	BioForm Stretch Soft Black Calf Panel, Double Layer of Neoprene, 4 + 2 Straps w/ Velcro Closure	X-Large	22	9
NDX-NB	BioForm Custom Stretch Soft Black Calf Panel, Layers of Neoprene, Straps w/ Velcro Closure	Custom	Template or Specs	

If color other than Black Neoprene is desired, replace the last two letters of the model number (NB) with those from Color Selection section.  
i.e.- NDS-NB becomes NDS-NE when color Eggplant is desired.

### B 3.7.1.3.- BioForm Non-Stretch Soft Black Calf Panels with Fastex Buckle Closure



Model	Description	Size	Dimension	
			W	H
VCS-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	Small	16	7
VCN-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	Medium	18	8
VCL-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	Large	20	8
VCO-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 2 Straps w/ Fastex Buckles	X-Large	22	9
VCX-V393	BioForm Custom Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), Straps w/ Fastex Buckles	Custom	Template or Specs	

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.-  
VCS- V393 becomes VCS-V322 when color Sun Yellow is desired.

***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

#### B 3.7.1.4.- BioForm Non-Stretch Soft Black Calf Panels with Velcro Closure (4+2 straps)



Model	Description	Size	Dimension	
			W	H
VDS-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	Small	16	7
VDM-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	Medium	18	8
VDL-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	Large	20	8
VDO-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), 4 + 2 Straps w/ Velcro Closure	X-Large	22	9
VDX-V393	BioForm Non-Stretch Soft Black Calf Panel, Covered w/ Vinyl (Naugahyde), Straps w/ Velcro Closure	Custom	Template or Specs	

If color other than Black Vinyl is desired, replace the last four digits of the model number (V393) with those from Color Selection section. i.e.- VDS- V393 becomes VDS-V322 when color Sun Yellow is desired.

#### B 3.7.2.- BioForm Solid Calf Panels

The BioForm Calf Panels are constructed from a ½” HDPE base, ½” Sunmate Foam, cover and mounting hardware. They are available with fixed mounting hardware as well as removable hardware.

##### B 3.7.2.1.- BioForm Solid Calf Panel, Fixed Mounting Hardware



Model	Description	Size	Dimension	
			W	H
CPS-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 13” & 14” Wheelchairs	Small	16	7
CPM-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 15” & 16” Wheelchairs	Medium	18	8
CPL-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 17” & 18” Wheelchairs	Large	20	8
CPO-NB	BioForm Solid Calf panel, ½” HDPE base, ½” Sunmate Foam covered in Black Neoprene, fixed hardware, for 19” & 20” Wheelchairs	X-Large	22	9
CPX-NB	BioForm Custom Solid Calf panel, fixed hardware	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section. i.e.- CPS-NB becomes CPS-NE when color Eggplant is desired.


***“If you don’t see it.... Contact Us”***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.7.2.2.- BioForm Solid Calf Panel, Removable Hardware



Model	Description	Size	Dimension	
			W	H
CRS-NB	BioForm Solid Calf panel, ½" HDPE base, ½" Sunmate Foam covered in Black Neoprene, removable hardware, for 13" & 14" Wheelchairs	Small	16	7
CRM-NB	BioForm Solid Calf panel, ½" HDPE base, ½" Sunmate Foam covered in Black Neoprene, removable hardware, for 15" & 16" Wheelchairs	Medium	18	8
CRL-NB	BioForm Solid Calf panel, ½" HDPE base, ½" Sunmate Foam covered in Black Neoprene, removable hardware, for 17" & 18" Wheelchairs	Large	20	8
CRO-NB	BioForm Solid Calf panel, ½" HDPE base, ½" Sunmate Foam covered in Black Neoprene, removable hardware, for 19" & 20" Wheelchairs	X-Large	22	9
CRX-NB	BioForm Custom Solid Calf panel, removable hardware	Custom	Template or Specs	

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.  
i.e.- CRS-NB becomes CRS-NE when color Eggplant is desired.

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345


© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.

### B 3.8.- BioForm Anterior Knee Blocks

BioDynamics' BioForm Anterior Knee Blocks are designed to help properly position the lower extremities while assisting in stabilizing the pelvis. They are available in 2 styles of mounting options; center mount (removable) and center and lateral mount (removable with lateral straps for added security).

Anterior Knee Blocks are not recommended for clients with hip or knee problems and should not be utilized to hold the pelvis to the rear of the wheelchair.

Available in small, medium and large. Custom sizes available by selecting the appropriate model #.



**Anterior Knee Blocks**

Model	Description	Size	Dimension			
			W	H	WK	DK
CKS-NB	BioForm Anterior Knee Block, Center Mount, 14" wide.	Small	14	4	3	1
CKM-NB	BioForm Anterior Knee Block, Center Mount, 16" wide.	Medium	16	4.5	4	1
CKL-NB	BioForm Anterior Knee Block, Center Mount, 18" wide.	Large	18	5	5	1.5
CLS-NB	BioForm Anterior Knee Block, Center & Lateral Mount, 14" wide.	Small	14	4	3	1
CLM-NB	BioForm Anterior Knee Block, Center & Lateral Mount, 16" wide.	Medium	16	4.5	4	1
CLL-NB	BioForm Anterior Knee Block, Center & Lateral Mount, 18" wide.	Large	18	5	5	1.5
CKX-NB	BioForm Custom Anterior Knee Block to Spec.,	Custom	Template or Specs			

If covering other than Black Neoprene is desired, replace the last two digits of the model number (NB) with those from Color Selection section.  
i.e.- CKS-NB becomes CKS-NE when color Eggplant is desired.

***"If you don't see it.... Contact Us"***

BioDynamics Ltd – 160 Terminal Dr. – Plainview, New York 11803 – [www.BioDynamics.us](http://www.BioDynamics.us)

E-Mail: [quote@BioDynamics.us](mailto:quote@BioDynamics.us) Fax: 516.349.9229 - 866.659.2345 Voice: 516.777.2222 - 866.669.2345

© BioDynamics Ltd.2010. All rights reserved. Total or partial reproduction prohibited, except if used as a Reference of BioDynamics Products.